

HyperLadder™ IV

Shipping: Room Temperature Catalog numbers
Exp. Date: See vial BIO-33029: 200 Lanes
Batch No.: See vial BIO-33030: 500 Lanes

Storage and stability:

HyperLadder IV is shipped at room temperature. Stored at room temperature it is stable for 6 months, stored at +4°C it is stable for 12 months and stored at -20°C it is stable for 24 months from date of purchase. Avoid repeated freeze thawing cycles.

Safety precautions:

Harmful if swallowed. Irritating to eyes, respiratory system and skin. Please refer to the material safety data sheet for further information.



A Meridian Life Science® Company

Features

- Easy size determination
- 10 bands from 100bp – 1013bp
- Ready-to-use format
- Optional mass determination

Description

HyperLadder IV is a ready-to-use molecular weight marker, especially designed for easy size determination. This ready-to-use format reduces handling steps and saves time; simply transfer HyperLadder IV from the vial to the gel.

HyperLadder IV produces a pattern of 10 regularly spaced bands, ranging from 100bp to 1013bp. To allow easy identification and orientation, the 1013bp band has the highest intensity.

Optionally the HyperLadder can be used for mass determination.

A 5x sample-loading buffer is supplied for your convenience. Under no circumstances should it be used to dilute/load ladder.

Components

Product	200 Lanes	500 Lanes
HyperLadder™ IV	2 x 500µl	5 x 500µl
5x Sample Loading Buffer	1ml	1ml

Protocol

Vortex the tube and use 5µl per lane.

Use 5x loading buffer for sample loading, under no circumstances should it be used for diluting the ladder. The blue loading buffer runs at approximately 4bp on a 1% gel.

Composition

This marker is composed of a restriction digest plus one or more PCR products.

Trouble Shooting

Ladders not sinking - vortex briefly.

Notes

1. This product insert is a declaration of analysis at the time of manufacture
2. HyperLadder is a trademark of Bioline
3. No licence is conveyed in respect of US 6,680,378
4. Research Use Only

Associated Products

Product	Pack size	Cat. No.
Agarose	100g	BIO-41026
Agarose Tablets	150g	BIO-41028
DNA Loading Buffer, Blue	2 x 1ml	BIO-37045
DNA Loading Buffer, Red	2 x 1ml	BIO-37068
Tri-Color DNA Loading Buffer	2 x 1ml	BIO-37070

Product Citations

1. Meredith, J., *et al.* *Plos One*. **6(1)**, doi:10.1371/journal.pone.0014587 (2011).
2. Dzahini-Obiatey, H. & Fox, R. *Afr. J. Biotechnol.* **9(5)**, 593-603 (2010).
3. Van den Broeke, A., *et al.* *BMC Genomics* 2010, **11**, 179 (2010).
4. Tivendale, K.A., *et al.* *Microbiol.* **155**, 450-460 (2009).
5. Garshasbi, M., *et al.* *Amer. J. Human Gen.* **82(5)**, 3783-3792 (2008).
6. Buonocore, F., *et al.* *Mol. Immunol.* **45(11)**, 3168-3177 (2008).
7. Beutin, L., *et al.* *J. Clin. Microbiol.* **43(4)**, 1552-1563 (2005).

Bioline Ltd
UNITED KINGDOM

Tel: +44(0)20 8830 5300
Fax: +44 (0)20 8452 2822

Bioline USA Inc.
USA

Tel: +1 508 880 8990
Fax: +1 508 880 8993

Bioline GmbH
GERMANY

Tel: +49(0)33 7168 1229
Fax: +49 (0)337168 1244

Bioline (Aust) Pty. Ltd
AUSTRALIA

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763