

CH3-Blue Competent Cells

Shipping: On Dry Ice Catalog numbers
Exp. Date: See vial BIO-85039 $\geq 10^8$ cfu/ μ g of pUC19
Batch No.: See vial BIO-85040 $\geq 10^9$ cfu/ μ g of pUC19

Storage and stability:

CH3-Blue Competent Cells are shipped on Dry Ice and can be stored for up to 6 months at -80°C .

Product Specifications:

Efficiency	Pack Size	Control DNA
$\geq 10^8$ cfu/ μ g of pUC19	1ml (10 x 100 μ l)	pUC19 (10pg/ μ l)
$\geq 10^9$ cfu/ μ g of pUC19	1ml (20 x 50 μ l)	pUC19 (10pg/ μ l)

Genotype:

F⁻ $\Delta mcrA \Delta(mrr-hsdRMS-mcrBC) \Phi 80lacZ \Delta M15 \Delta lacX74 recA1 endA1 ara \Delta 139 \Delta(ara, leu)7697 galU galK \lambda-rpsL (Str^R)$

Notes:

1. This product insert is a declaration of analysis at the time of manufacture.
2. Research Use Only.

Store at -80°C



DATA SHEET

Features

- Lacks *mcrA*, *mcrBC*, *mrr* and *hsdRMS* restriction systems
- Available in two efficiencies: $\geq 10^8$ or $\geq 10^9$ cfu/ μ g of pUC19

Applications

- Cloning of methylated DNA
- Ideal for subcloning and generating cDNA libraries
- Blue/white color screening

Description

CH3-Blue Chemically Competent Cells are a highly efficient derivative of *E. coli* K12, ideal for the construction of cDNA libraries using plasmid derived vectors. To facilitate the cloning of DNA that contains methylcytosine or 5-hydroxymethylcytosine, CH3-Blue lacks the *E. coli* restriction systems *mcrA*, *mcrBC*, *mrr* and *hsdRMS*. The *lacZ* mutation allows blue/white color screening and α -complementation of recombinants. The *recA1* and *endA1* markers minimize recombination events and improve the quality and yield of plasmid DNA.

Suggested Transformation Procedure for Optimal Results:

1. Remove cells from -80°C and let thaw on wet ice.
2. Gently mix cells by lightly flicking tube. Aliquot ~ 50 - $100\mu\text{l}$ of cells into chilled, 17 x 100mm polypropylene tube(s), e.g., Falcon 2059. Unused cells may be refrozen, but a small drop in efficiency may result. For optimal recovery, refreeze cells in a dry ice/ ethanol bath prior to -80°C storage.
3. Add DNA solution ($\leq 5\mu\text{l}$ per $50\mu\text{l}$ cells) to cell suspension and gently swirl tube(s) for a few seconds to mix. If a control is desired, repeat this step with $2\mu\text{l}$ of the provided pUC19 in a separate tube.
4. Incubate on ice for 30 minutes.
5. Place tube(s) in 42°C water bath for ~ 30 to 45 seconds without shaking. For $50\mu\text{l}$ aliquots in Falcon 2059 tubes, 30 seconds is recommended for maximum efficiency.
6. Replace tube(s) on ice for ~ 2 minutes.
7. Dilute transformation reaction(s) to 1ml by addition of 900-950 μl SOC. SOC Medium: 2% Tryptone, 0.5% Yeast Extract, 0.4% glucose, 10mM NaCl, 2.5mM KCl, 10mM MgCl_2 & 10mM MgSO_4 .
8. Shake tube(s) ~ 200 rpm for 60 minutes at 37°C .
9. Plate by spreading 5- $200\mu\text{l}$ of cell transformation mixture on LB agar plates containing appropriate antibiotic and incubate overnight at 37°C .

When performing the pUC19 control transformation, plate $5\mu\text{l}$ of the transformation mixture on a LB agar plate containing 100 $\mu\text{g}/\text{ml}$ ampicillin. To facilitate cell spreading, place a pool of SOC (100 μl) onto surface of plate prior to addition of transformation mixture.

Transformation Efficiency Calculation for Control DNA

$$\text{Transformation Efficiency (cfu}/\mu\text{g pUC19 DNA)} = \frac{\text{\# colonies (colony forming units)}}{\text{pg pUC19 transformed}} \times \frac{10^6 \text{ pg}}{\mu\text{g}} \times \frac{\text{Final volume } (\mu\text{l}) \text{ of transformation mix}}{\text{Volume plated } (\mu\text{l})}$$

For example:

If 40 colonies were obtained after transforming 20pg of pUC19 and plating $5\mu\text{l}$ of the final 1ml transformation mixture, the calculated transformation efficiency would be:

$$\frac{40 \text{ cfu}}{20 \text{ pg pUC19}} \times \frac{10^6 \text{ pg}}{\mu\text{g}} \times \frac{1000 \mu\text{l}}{5 \mu\text{l}} = 4 \times 10^8 \text{ cfu}/\mu\text{g pUC19}$$

Associated Products:

Product Name	Pack Size	Cat No
T4 DNA Ligase	500 Units	BIO-27026
Quick-Stick Ligase	50 Reactions	BIO-27027
IPTG	5g	BIO-37036
X-GAL	1g	BIO-37035

Product Citations:

1. Thompson, K. M. *et al. FEMS Micro. Lett.* **305(2)**, 143-7 (2010).