

INSTRUCTIONS FOR USE

**Babesia canis IFA
IgG Antibody Kit**

Catalog Number: CBG-100

Size: 100 test

Storage: 2-8°C

An Indirect fluorescence immunoassay for the detection of IgG class antibody against *Babesia canis* in canine serum or plasma

For in-vitro diagnostic use only



INTENDED USE

The *Babesia canis* IFA IgG Antibody kit is intended for the detection and semi-quantitation of IgG class canine antibody to *Babesia canis*. This kit is designed for use as an aid in the diagnosis of canine infection by this pathogen.

SUMMARY AND EXPLANATION OF TEST

The IFA slides in this kit contain canine erythrocytes, 20-30% of which are infected with *Babesia canis*. Test sera are diluted to screening dilution in phosphate-buffered saline (PBS) and incubated in the individual slide wells to allow reaction of serum antibody with the solid-phase antigens. The slides are then washed to remove unreacted serum proteins, and an DyLight 488-labeled anti-canine IgG (Conjugate) is added to label the antigen-antibody complexes. After further incubation, the slides are washed again to remove unreacted Conjugate. The resulting reactions are visualized using standard fluorescence microscopy, where a positive reaction is seen as sharply defined apple-green fluorescent *Babesia* inclusions (merozoites) within the erythrocytes. A negative reaction is seen as either red-counterstained cells or fluorescence unlike that seen in the Positive Control well. Positive reactions may then be retested at higher dilutions to determine the highest reactive or endpoint dilution.

GENTHAUR

REAGENTS

IFA Ag x 12

Substrate Slides (10)

10 x 10-well masked slides containing canine erythrocytes infected with *Babesia canis*. Slides are fixed (inactivated) and packaged under vacuum.

CONJ FITC

Conjugate, 2.5 mL

Yellow cap dropper bottle contains affinity-purified DyLight 488-labeled rabbit anti-canine IgG (heavy chain) with bovine serum albumin and Evans' blue counterstain.

CONT +

Positive Control, 0.5 mL

Blue cap dropper bottle contains reactive canine serum at a 1:50 screening dilution. Endpoint titer is 1:400

CONT -

Negative Control, 0.5 mL

Red cap dropper bottle contains non-reactive canine serum at a 1:50 screening dilution

MM

Mounting Medium, 1 mL

White cap dropper bottle contains glycerol (50% v/v) in PBS

BUF WASH PBS

PBS, 1 liter

Add supplied powder to 1 liter purified water to produce PBS.

Warnings

Since no testing can assure the absence of infectious agents these reagents, as well as all serum specimens and equipment coming in contact with these specimens, should be handled with good laboratory practices to avoid skin contact and ingestion. The substrate slides are prepared with chemically inactivated antigens. However, the slides should be considered potentially infectious and handled accordingly.

Storage

Kit components should be stored at 2-8°C or colder. Bring them to room temperature (20-25°C) before opening bottles or slide envelopes.

SPECIMEN COLLECTION

Allow blood sample to clot and separate sera by centrifugation. Transfer sera aseptically to tightly closing sterile containers. Store at 2-8°C. If testing is to be delayed longer than 5 days, storage at -20°C or colder is recommended. Acute specimens should be drawn at the onset of illness; convalescent specimens should be obtained at two and four week intervals to check for titer changes.

PROCEDURE

This kit supplies sufficient reagents and materials for 100 determinations.

Materials Required But Not Supplied

- Purified (distilled or deionized) water
- Clean 250 or 500 mL wash bottle for PBS
- 12x75 mm test tubes or microtiter plate for preparing serum dilutions
- Precision pipette(s) in microliter range for making and delivering serum dilutions
- 24 x 50 mm glass cover slips
- Fluorescence microscope with filter system for FITC (maximum excitation wavelength 490 nm, mean emission wavelength 530 nm) and 400X magnification
- 37° water bath or incubator
- Humidity chamber for slide incubation steps

Precautions

- Do not use components past expiration date.
- Conjugate is photosensitive and is packaged in opaque plastic for protection. Store in the dark at 2-8°C and return to storage immediately after use.
- Conjugate contains Evans' Blue dye, which may be carcinogenic. Avoid contact with skin.

ASSAY PROCEDURE

Before opening slide pouches, allow them to attain ambient temperature.

7. Wash slide as in step 5, above. Then add 2-3 drops Mounting Medium to each slide and cover slip.
9. Read the stained substrate slide at 400X magnification. Slide may be stored at 2-8°C in the dark for up to 24 hours.

QUALITY CONTROL

The Negative Control and dilutions of the Positive Control should be assayed with each daily run. The Negative Control well is an example of a non-reactive serum, with either uniform red counterstain or slight (less than 1+), but uniform, greenish staining. The Positive Control wells should give an endpoint titer from 1:100 to 1:800. The fluorescence intensity at 1:400 may be used as the cut-off level required for a test reaction to be called positive. If either of the Controls does not react as specified, the assay run should be considered void. Reagent components and procedural steps should be rechecked, and the assay repeated from the beginning.

The Negative Control well is an example of fluorescence patterns that are to be considered negative. If bright staining is seen in this well, similar to that seen in the Positive Control wells, there has been a breakdown in technique and the assay must be repeated.

INTERPRETATION OF RESULTS

A positive reaction appears as brightly fluorescent (at least 1+) sharp, regular stained inclusions within the erythrocytes, while the Negative Control well is an example of fluorescence patterns that are to be considered negative. The Positive Control wells should give an endpoint titer of 1:200 to 1:800. The fluorescence intensity at 1:400, however, may be used as the cut-off level required for a patient reaction to be called positive.

The size, appearance and density of the inclusions must be compared with the Positive and Negative Control reactions. Patterns of reactivity different than those seen in the Positive Control must be considered non-specific.

PATIENT SPECIMENS

GENT

1. Prepare 1:50 screening dilutions (1 part test serum with 49 parts PBS) for all untested serum specimens. For sera found positive on a previous assay run, prepare serial two-fold dilutions in PBS, starting with 1:50. Acute-convalescent pairs should be compared by assaying all dilutions in parallel.
2. Prepare dilutions of the Positive Control in PBS to include one dilution above the stated endpoint and one dilution below the stated endpoint (1:200-1:800). Note: These controls are bottled at screening dilution (1:50), so dilute no more than 8-fold further to reach endpoint.
3. For each serum or Control dilution to be tested, add 10 μ L to one slide well. For each assay, include the Negative Control, Positive Control and dilutions of the Positive Control (step 2).
4. Place slides into a humidity chamber and incubate in water bath or incubator for 30 minutes at $37^{\circ}\pm 0.5^{\circ}\text{C}$.
5. Rinse slide wells with gentle stream of PBS from the wash bottle three (3) times, shaking PBS from the slide into a sink between each wash. Go directly to the next step without allowing slide wells to dry.
6. To each slide well, add 10 μ L Conjugate then return slide to the humidity chamber for 30 minutes incubation in the water bath or incubator at $37^{\circ}\pm 0.5^{\circ}\text{C}$. Incubation should be in the dark to protect the photosensitive conjugate.

Positive at 1:50: Single IgG titers of 1:50 and greater are considered to reflect recent or active infection.

Negative at 1:50: Report as negative for antibody.

EXPECTED VALUES

The prevalence of *Babesia canis* antibodies varies depending upon the geographic region and population being tested.

Revised 11/2006

GENZYME