

Store at -20°C

Cat. No.	Conc.	Quantity
G078	5U/μl	500U

Description

Precision™ Taq DNA Polymerase is a thermostable enzyme that replicates DNA at 70°C to 80°C. It catalyzes the polymerization of nucleotides into duplex DNA in the 5' to 3' direction in the presence of magnesium. The enzyme has a molecular weight of approximately 90,000 daltons. Unlike Taq DNA Polymerase, Precision™ Taq DNA Polymerase exhibits 3' to 5' exonuclease (proofreading) activity, that enables the polymerase to correct nucleotide incorporation errors. The error rate of Precision™ Taq DNA Polymerase in PCR is 2.6x10⁻⁶ per nucleotide per cycle. The enzyme can be used in PCR applications that demand high fidelity and can amplify templates up to 3-5kb in length. The amplified PCR products are blunt-ended.

Source

E.coli cells carrying a cloned *pol* gene from *Pyrococcus Furiosus*.

Components

Precision™ DNA Polymerase	100μl
10X PCR Buffer (-MgCl ₂)	1.0ml
25mM MgCl ₂	1.0ml

Enzyme Storage Buffer

50mM Tris-HCl (pH 8.2 at 25°C), 0.1mM EDTA, 1mM DTT, 0.05% CHAPS and 50% glycerol.

10X PCR Buffer

200mM Tris-HCl (pH 8.8), 100mM KCl, 100mM (NH₄)₂SO₄, 1% TritonX-100 and enzyme stabilizer. 25mM MgCl₂ supplied separately.

Unit Definition

One unit of the enzyme catalyzes the incorporation of 10 nanomoles of deoxyribo-nucleotides into a polynucleotide fraction (adsorbed on DE-81) in 30 mins at 72°C.

Basic PCR Protocol

The following basic protocol serves as a general guideline and a starting point for any PCR amplification. Optimal reaction conditions (incubation times and temperatures, concentration of Taq Polymerase, primers, MgCl₂ and template DNA) vary and need to be optimized.

PCR reactions should be assembled in a DNA-free environment. DNA sample preparation, reaction mixture assemblage and the PCR process, in addition to the subsequent reaction analysis, should be performed in separate areas.

A control reaction, omitting template DNA, should always be performed to confirm the absence of contamination.

1. Add the following components to a sterile 0.2ml PCR tube sitting on ice.

Components	Volume
DNA Template	~100ng
Forward Primer (10μM)	1μl
Reverse Primer (10μM)	1μl
10X PCR Buffer	5μl
25mM MgCl ₂	3μl
dNTP Mixture (10mM each)	1ul
Precision™ Taq (5U/μl)	0.5-1μl
ddH ₂ O	up to 50μl

We recommend preparing a mastermix for multiple reactions to minimize reagent loss and enable accurate pipetting.

2. Mix contents of tube and centrifuge briefly.
3. Incubate tube in a thermal cycler at 94°C for 3 mins to completely denature the template.
4. Perform 30-40 cycles of PCR amplification as follows:

Denature: 94°C for 30 sec
Anneal: 55°C for 30 sec
Extend: 72°C for 1 min/1kb template

5. Incubate for an additional 5 mins at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.
6. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide or SafeView™ staining. Use appropriate molecular weight standards.

*This product is distributed for laboratory research only.
CAUTION: Not for clinical use.
For technical questions about this product, phone the ABM helpline at 1-866-571-7226 or visit our website at www.abmGood.com.*

CERTIFICATE OF ANALYSIS