





Total RNA
Extraction System



Plant Total RNA
Extraction System



Viral RNA
Extraction System



Ultrapure Plasmid
Extraction System


Mini-M™
Plasmid DNA
Extraction System


Gel-M™
Gel Extraction
System


PCR-M™
Clean Up System


Genomic DNA
Extraction System


Plant
Genomic DNA
Extraction System

Total RNA Extraction Miniprep System

User's Guide



Service

Viogene regards it very important to provide satisfactory service to our every customer. In order to guarantee the best quality of our products, we value our customers' comments and suggestions on our services, or the performance, new applications, and techniques of our products. If there is any question or comment concerning the use of our products, please do not hesitate to contact our Technical Service Department by phone, e-mail, or fax, or to contact your local sales representatives. Our experienced staffs and researchers are pleased to provide you with technical help and advice. If you have problems on attaining the expected performance with our products, please contact our Technical Service Department for technical advice. If any product fails to perform properly not due to incorrect handling, please contact us or your local sales representatives for assistance.

Quality Control

We strictly require good quality control of our products by regular testing of each lot to maintain a satisfactory yield of DNA or RNA. Testing results of all lots of each product are documented. Any inquiry to access them is welcome.

Table of Contents

Contents	Page no.
Description	5
Downstream Applications	5
Product Contents	6
Shipping and Storage	6
Important Notes	7-8
Protocol	
I. Animal Tissue Protocol	9-10
II. Animal Cell Protocol	10-11
III. Animal Cell Cytoplasm Protocol	11
IV. Bacteria Protocol	12
V. Removal of genomic DNA in eluted total RNA	13
VI. Application of Shearing Tube	13
Troubleshooting Guide	14-16
Viogene's Hints	17-18
Viogene Products	19



Description

Viogene Total RNA Extraction Miniprep System provides an economical method to purify total RNA from various samples such as cultured cells, tissues, and bacteria. A simple silica-membrane spin-column method can isolate total RNA without need of performing time-consuming phenol/chloroform extraction and ethanol precipitation. Total RNA longer than 200 nucleotides are isolated, while small RNA such as 5.8S RNA, 5S RNA, and tRNA, which make up 15-20% of the total RNA, are excluded.

Sample Preparation Guide (Preparation Time: ~30 minutes)

Sample	Recommended amount of sample used	Maximum Yield (μg)	
Animal cells	NIH-3T3	1×10^6 cells	12
	HeLa	1×10^6 cells	15
	COS-7	1×10^6 cells	30
	LMH	1×10^6 cells	12
Animal tissues (mouse/rat)	Embryo	10 mg	30
	Heart	10 mg	10
	Brain	10 mg	10
	Kidney	10 mg	35
	Liver	10 mg	45
	Spleen	10 mg	35
	Lung	10 mg	10
Bacteria	Thymus	10 mg	45
	<i>E. coli</i>	1×10^9 cells	65
	<i>B. subtilis</i>	1×10^9 cells	40

Downstream Applications

- * Northern blotting
- * RT-PCR
- * Poly A⁺ RNA selection
- * cDNA synthesis
- * Primer extension
- * *in vitro* translation

Product Contents

	GR1001 (50 preps) (Cat. No.)	GR1002 (250 preps) (Cat. No.)
RX Buffer	36 ml (GR1001RX01)	200 ml (GR1002RX01)
WF Buffer	30 ml (KR1001SWF)	150 ml (KR1002SWF)
WS Buffer	15 ml* (KR1001SWS)	45 ml** (KR1002SWS)
RNase-free ddH ₂ O	1.5 ml x 2	15 ml
Total RNA Mini Column (RNase-free)	50 pieces	250 pieces
Collection Tube (RNase-free)	50 pieces	250 pieces
1.5-ml Elution Tube (RNase-free)	50 pieces	250 pieces
Protocol	1	1

* For GR1001 (50 preps), add 60 ml of 98-100% ethanol into WS Buffer bottle when first open.

** For GR1002 (250 preps), add 180 ml of 98-100% ethanol into WS Buffer bottle when first open.

Buffers are available for separate purchase. Please refer to the Cat. No. listed above for ordering.

Shearing Tube
(RNase-free) 10 pieces/pk
(GRS1001)

Shearing Tube is not provided in this system. Separate order is required.

Shipping and Storage

All components of VioGene Total RNA Extraction Miniprep System are stable at room temperature (20-25°C) for one year.

Important Notes

Please read the following notes before starting the procedures.

1. All plasticware and containers should be treated properly to make sure RNase-free. Gloves should be worn when handling RNA.
2. Buffers contain irritants. Appropriate safety apparels such as gloves and lab coat should be worn to protect from skin contact.
3. All procedure should be done at room temperature (20-25°C).
4. Pipet a required volume of RX Buffer into another tube and add 10 μ l β -mercaptoethanol (β -ME) per 1 ml RX Buffer before use.
5. For GR1001 (50 preps), add 60 ml of 98-100% ethanol into WS Buffer bottle **when first open**. For GR1002 (250 preps), add 180 ml of 98-100% ethanol into WS Buffer bottle **when first open**. Ethanol is provided by the user.
6. Do **not** use more than the suggested maximum amount of sample (refer to **VioGene's Hints**, No. 1, page 17).

- Complete disruption and homogenization of sample is essential for total RNA extraction.
- After each vortexing step, when the tube is opened, **briefly centrifuge** the tube to bring down the sample attached inside the cap to avoid generation of aerosols and contact with sample.
- When sample or buffer is added into the column, avoid touching the rim. This is to prevent cross contamination of samples when handling the columns.
- All centrifugation steps except cell pelleting are done at full speed (10,000 x g or 13,000-14,000 rpm) in a microcentrifuge.
- Some genomic DNA (and plasmid DNA, if any) will also be copurified with RNA. DNase treatment is therefore required when DNA-free RNA is desired. DNase can then be removed by phenol/chloroform extraction (refer to **Protocol** for "Removal of genomic DNA in eluted total RNA", page 13).

Protocol

Please refer to the Table of Contents on page 4 to choose the appropriate protocol according to the kind of sample used.

I. Animal Tissue Protocol

- Add 350 μ l RX Buffer (β -ME added) to 10 mg of liquid-nitrogen-frozen or fresh tissue. Disrupt and homogenize the sample by grinding and shearing using 20-G needle syringe or Viogene's Shearing Tube.
- Centrifuge the lysate for 5 minutes to spin down insoluble materials and use only the supernatant in the following steps.
- Determine the final volume of the supernatant. Add an equal volume of 70% ethanol to the clear lysate and mix by vortexing.
- Place a Total RNA Mini Column onto a Collection Tube. Add 700 μ l of the ethanol-added sample (including any precipitate) into the column. Centrifuge for 30-60 seconds. Discard the flow-through.
- Wash the column once with 0.5 ml WF Buffer by centrifuging for 30-60 seconds. Discard the flow-through.

Add 10 μ l β -mercaptoethanol (β -ME) per 1 ml RX Buffer.

If use 20 mg tissue, add 700 μ l RX Buffer to ensure complete sample lysis.

If Viogene's Shearing Tube is used, refer to page 13 for "Application of Shearing Tube".

A gelatinous layer of substance may form at the bottom of the tube after centrifugation. Avoid taking it up when take out the supernatant. If this layer is not well pelleted, further centrifuge the tube for a few minutes more.

70% ethanol should be prepared using DEPC-treated ddH₂O.

Repeat this step for the rest of the sample.

If some sample still retains in the column, repeat centrifugation until all sample pass the column.

6. Wash the column once with 0.7 ml WS Buffer by centrifuging for 30-60 seconds. Discard the flow-through.

Make sure that ethanol has been added into the WS Buffer bottle when first open.

7. Centrifuge the column for another 3 minutes to remove ethanol residue.

Residual ethanol can affect the quality of RNA and inhibit subsequent enzymatic reactions such as reverse transcriptase reaction. If necessary, centrifuge the column for a few minutes more to remove all the ethanol.

8. Place the column onto a 1.5-ml RNase-free Elution Tube. Add 30-50 μ l RNase-free ddH₂O (provided) onto the membrane.

For effective elution, make sure that the elution solution is dispersed onto the center of the membrane.

Eluting the column twice can result in a higher RNA recovery (refer to **Viogene's Hints**, No. 4, page 16).

9. Stand the column for 1 minute, and centrifuge for 1-2 minutes to elute total RNA.

10. Store RNA at -70°C.

I. Animal Cell Protocol

1. Pellet 1 to 5 x 10⁶ cells by centrifuging at 300 x g for 5 minutes. Remove all the supernatant.

Any residual supernatant present will affect cell lysis by RX Buffer.

2. Disrupt cells by adding 350 μ l RX Buffer (β -ME added) to the cell pellet and vortex the sample. Homogenize the sample by using 20-G needle syringe or Viogene's Shearing Tube.

Add 10 μ l β -mercaptoethanol (β -ME) per 1 ml RX Buffer.

Add 700 μ l RX Buffer when use 1 x 10⁷ cells.

If Viogene's Shearing Tube is used, refer to page 13 for "Application of Shearing Tube".

3. Follow the **Animal Tissue Protocol** starting from Step 2 on page 9.

III. Animal Cell Cytoplasm Protocol

1. Prepare cytoplasm lysate:

Prepare cell lysis buffer (20 mM Tris-HCl, pH 8.0; 1 mM MgCl₂; 0.5% NP-40). Keep at 4°C.

Only **fresh** cells are used for preparing cytoplasm lysate.

Total RNA extracted from cytoplasm lysate are of minimum genomic DNA contamination.

a. Pellet 5 x 10⁶ to 1 x 10⁷ fresh cells by centrifuging at 300 x g for 5 minutes. Remove all the supernatant.

b. Add 180 μ l cell lysis buffer (4°C) to the cell pellet, resuspend and lyse cells by **gentle** pipetting. Incubate the lysate on ice for 5 minutes.

c. Centrifuge the lysate at 300 x g at 4°C for 3 minutes, transfer the supernatant to a new tube and discard the pellet. Use the supernatant (lysate) in the following steps.

2. Add 600 μ l RX Buffer (β -ME added) to the lysate and mix by vortexing.

Add 10 μ l β -mercaptoethanol (β -ME) per 1 ml RX Buffer.

3. Add 450 μ l 98-100% ethanol to the sample and mix by vortexing.

4. Follow the **Animal Tissue Protocol** starting from Step 4 on page 9.

IV. Bacteria Protocol

1. Pellet up to 1×10^9 bacterial cells by centrifuging at $5,000 \times g$ (7,500 rpm) for 5 minutes. Remove all the supernatant.
2. Resuspend cells in 100 μ l TE buffer by vortexing.
3. Add lysozyme to a final concentration of 500 μ g/ml for **Gram-negative bacteria**; 2 mg/ml for **Gram-positive bacteria**, and incubate at room temperature for 5 to 10 minutes to digest the cell wall.
4. Add 350 μ l RX Buffer to the sample and mix by vortexing.
5. Centrifuge lysate for 5 minutes to spin down insoluble materials and use only the supernatant in the following steps.
6. Add 250 μ l 98-100% ethanol to the sample and mix by vortexing.
7. Follow the **Animal Tissue Protocol** starting from Step 4 on page 9.

Lysozyme is provided by user.

Add 10 μ l β -mercaptoethanol (β -ME) per 1 ml RX Buffer.

A gelatinous layer of substance may form at the bottom of the tube after centrifugation. Avoid taking it up when take out the supernatant. If this layer is not well pelleted, further centrifuge the tube for a few minutes more.

V. Removal of genomic DNA in eluted total RNA

1. Incubate total RNA with RNase-free DNase I (1 unit DNase I per μ g RNA) in 50 mM Tris-HCl (pH 7.5), 10 mM $MnCl_2$, and 50 μ g/ml BSA at 37°C for 15-30 minutes.
2. Remove DNase I by adding an equal volume of phenol:chloroform (1:1) and mix well. Centrifuge for 5 minutes. Transfer the upper aqueous layer to a new eppendorf tube.
3. Add $1/10$ volume of 3 M sodium acetate (pH 5.2) and 1 volume of ice-cold isopropanol to the solution and mix well. Chill on ice for 30 minutes.
4. Centrifuge for 10 minutes at 4°C. Wash RNA pellet twice with 1 ml of 70% ethanol and recentrifuge.
5. Remove all supernatant. Air dry RNA pellet. Redissolve RNA in RNase-free ddH_2O .

VI. Application of Shearing Tube

Shearing Tube is designed for simple and fast homogenization of tissue and cell lysate. The lysate is loaded into a Shearing Tube sitting in a 2-ml Collection Tube and centrifuge the tube for 1-2 minutes at full speed ($10,000 \times g$ or 13,000-14,000 rpm) in a microcentrifuge. When collecting homogenized lysate from the Collection Tube, avoid pipetting any debris and pellet formed at the bottom of the tube.

Troubleshooting Guide

Problem	Possible Reason	Solution
Column is clogged when passing the sample	Sample lysate contains insoluble residues or/and gelatinous substance	After sample lysis, centrifuge the sample at full speed for 5 minutes or more and only use the supernatant.
	Sample lysate is very viscous because too much sample is used	Reduce the sample amount or increase the volume of RX Buffer and ethanol proportionally.
Little or no RNA eluted	Insufficient disruption or homogenization	Reduce the amount of starting sample and perform more disruption and homogenization to the sample.
	Column is clogged	Reduce the amount of starting sample and perform more disruption and homogenization. Centrifuge the lysate to remove insoluble materials and use only the supernatant.
	RNA is not completely eluted because RNase-free ddH ₂ O does not penetrate into the membrane	Add ddH ₂ O onto the center of the membrane and stand the column for 5 minutes. If ddH ₂ O still retains on the membrane, pulse centrifuge the column for a few seconds to drag ddH ₂ O into the membrane.
	RNA is degraded	Flash freeze fresh samples in liquid nitrogen and store at -80°C if not used immediately. Improper handling (such as thawing) of the sample or storing the sample at -20°C will cause RNA degradation.

Problem	Possible Reason	Solution
Little or no RNA eluted	RNase contamination	Treat bench surface before use. Use RNase-free solutions, plasticware and glassware.
	No ethanol or ethanol of incorrect amount is added to the sample lysate	Determine the final volume of sample lysate obtained. Add ethanol of correct volume and concentration as indicated in the protocol.
DNA contamination	DNA is copurified with RNA	Use RNase-free DNase to treat the eluted RNA sample. DNase can then be removed by phenol/chloroform extraction (refer to Protocol , page 13).
		Minimize DNA copurification by extracting total RNA from cytoplasm lysate prepared from fresh cultured cells .
A ₂₆₀ /A ₂₈₀ ratio of eluted total RNA is low	Use ddH ₂ O of acidic pH to dilute RNA samples for spectrophotometric analysis	Use 10 mM Tris-HCl of pH 7.5 or TE buffer to dilute the RNA samples (refer to Viogene's Hints , No. 3, page 18).
	Proteins in the sample are not completely denatured	Too much sample is used. Reduce the sample amount or increase the volume of RX Buffer and ethanol proportionally.

Troubleshooting Guide

Problem	Possible Reason	Solution
A_{260}/A_{280} ratio of eluted total RNA is low	DNA is copurified with RNA Eluted RNA carries contaminants	Refer to Solution section of Problem - "DNA contamination". Wash the column twice with 0.7 ml WS Buffer.
RNA appears smearing and degraded	Sample is stored and handled improperly	Flash freeze fresh samples in liquid nitrogen and store at -80°C if not used immediately. Improper handling (such as thawing) of the sample or storing the sample at -20°C will cause RNA degradation.
	Cell samples were harvested from an old or over-grown culture	RNA of good quality is only expected from a healthy cell or bacterial culture.
	RNase contamination	Treat bench surface before use. Use RNase-free solutions, plasticware and glassware.
	Gel electrophoresis is performed in running buffer or tank contaminated with RNase	Use fresh running buffer prepared from DEPC-treated ddH ₂ O and properly-cleaned tank for electrophoresis.
Poor performance in downstream applications	Eluted RNA carries ethanol residue	After wash with WS Buffer, do discard the flow-through, and centrifuge the column for another 3 minutes. If necessary, centrifuge the column for a few minutes more to remove all the ethanol.

Viogene's Hints

1. Low yield or purity of total RNA is usually due to incomplete lysis of the sample. Starting with a maximum amount of samples does NOT usually give the best yield of RNA. Instead, it always results in incomplete sample lysis and degradation of proteins as well as reduces RNase inactivation effectiveness of RX Buffer, thus affecting mRNA intactness and making extraction of all the total RNA from the sample unfeasible. Besides, it also requires subsequent removal of insoluble residues and yields viscous sample lysate. When the lysate is too viscous, it not only has difficulty in passing the column, but also indicates the presence of an abundant amount of contaminants such as proteins and salts. Contaminants of high amount not only affect RNA binding, but also may not be washed off completely, leading to carry over to the eluted RNA. Therefore, a good quality and yield of RNA is only expected when a sample is **completely** lysed. We advise starting with half of the maximum amount of sample suggested. When there is no problem in completely lysing the sample and passing the lysate through the column, amount of the sample to be applied can be increased gradually in the subsequent preparations.
2. Though Viogene Total RNA Extraction Miniprep System can also be used to extract total RNA from plant samples, Viogene Plant Total RNA Extraction Miniprep or Maxiprep System is recommended for plant samples for better sample lysis effect.

3. Using ddH₂O of acidic pH (5.0-6.0) to dilute DNA and RNA samples for spectrophotometric analysis will significantly decrease A₂₆₀/A₂₈₀ ratio of the sample (Wilfinger et al., 1997). 10 mM Tris-HCl of pH 7.5 or TE buffer should be used to dilute the samples.
4. When a more concentrated RNA sample is desired, elute with 30 µl RNase-free ddH₂O. Eluting the column **twice** with 30-50 µl RNase-free ddH₂O can result in a higher RNA recovery, especially when expected RNA yield is more than 30 µg. The two eluants can be collected into the same or separate tube.

Reference: Wilfinger, W. W., Mackey, K., and Chomczynski, P. 1997. Effect of pH and ionic strength on the spectrophotometric assessment of nucleic acid purity. *Biotechniques* **22**:474-481.

VioGene Products

Product	Cat. No.	Package Size	Sample	Expected Yield
Mini-M Plasmid	GF1001	50	1-5 ml culture	up to 20 µg
	GF1002	250		
Midi-V100 Plasmid	GDV1001	25	25-100 ml culture	up to 100 µg
	GDV1002	50		
Maxi-V500 Plasmid	GMV1001	10	100-250 ml culture	up to 500 µg
	GMV1002	25		
Blood & Tissue Genomic DNA Mini	GG1001	50	200 µl whole blood	up to 10 µg
	GG1002	250		
Blood Genomic DNA Midi	GGD1001	20	1 ml whole blood	up to 50 µg
	GGD1002	100		
Blood Genomic DNA Maxi	GGM1001	10	5 ml whole blood	up to 300 µg
	GGM1002	50		
Plant Genomic DNA Mini	GPG1001	50	100 mg tissue	up to 40 µg
	GPG1002	250		
Plant Genomic DNA Maxi	GPGM1001	20	1 g tissue	up to 1 mg
Total RNA Mini	GR1001	50	10-20 mg tissue	10-45 µg
	GR1002	250		
Total RNA Midi	GRD1001	10	0.1-0.2 g tissue	200-450 µg
	GRD1002	50		
Total RNA Maxi	GRM1001	6	0.5-1 g tissue	1-5 mg
	GRM1002	24		
Viral RNA Mini	GVR1001	50	150 µl body fluid	up to 90% recovery
	GVR1002	250		
Plant Total RNA Mini	GPR1001	50	100 mg tissue	up to 100 µg
	GPR1002	250		
Plant Total RNA Maxi	GPRM1001	10	1 g tissue	up to 1 mg
Gel-M Gel Extraction	EG1001	50	50-200 mg agarose gel	50-80% recovery (100 bp-10 kb)
	EG1002	250		
PCR-M Clean Up	PF1001	50	10-100 µl DNA	up to 95% recovery (100 bp-10 kb)
	PF1002	250		
VioTag DNA Polymerase	VT1001	500 U (5 U/µl)	10X PCR Buffer containing 20 mM MgCl ₂	
VioTwinPack Kit	VTP1001	500 U VioTag DNA Polymerase (5 U/µl)	10X PCR Buffer containing 20 mM MgCl ₂	40 mM of dNTP mix (10 mM each)
Clear-band Agarose	AG0050	50 g		
	AG0100	100 g		