

Hydrogen Peroxide Assay Kit (Catalog #K265-200; 200 reactions; Store kit at -20°C)

I. Introduction:

Hydrogen Peroxide is a reactive oxygen metabolic byproduct that serves as a key regulator for a number of oxidative stress-related states. Functioning through NF-KB and other factors, hydroperoxide-mediated pathways have been linked to asthma, inflammatory arthritis, atherosclerosis, diabetic vasculopathy, osteoporosis, neurodegenerative diseases, Down's syndrome and immune system diseases. BioVision's Hydrogen Peroxide Assay Kit is a highly sensitive, simple, direct and HTS-ready colorimetric and fluorometric assay for measuring H₂O₂ in biological samples. In the presence of Horse Radish Peroxidase (HRP), the OxiRed Probe reacts with H₂O₂ to produce product with color ($\lambda_{max} = 570 \text{ nm}$) and red-fluorescent (Ex/Em=535/587 nm). The kit can perform 200 reactions by fluorometric method or 100 reactions by colorimetric method. The detection limit can be as low as 2 pmol per assay (or 40 nM concentration) of H₂O₂ in the sensitive fluorometric assay.

II. Kit Contents:

Components	K265-200	Cap Code	Part No.
H ₂ O ₂ Assay Buffer	25 ml	WM	K265-200-1
OxiRed™ Probe	1 vial	Red	K265-200-2
Dimethylsulfoxide (DMSO, anhydrous)	0.4 ml	Brown	K265-200-3
HRP	1 vial	Green	K265-200-4
H ₂ O ₂ Standard (0.88M)	0.1 ml	Yellow	K265-200-5

III. Storage and Handling:

Warm the assay buffer to room temperature before use. Briefly centrifuge vials prior to opening. Read the entire protocol before performing the assay.

IV. Reagent Reconstitution and General Consideration:

OxiRed™ Probe: Dissolve in 220 μl DMSO (provided), pipeting up and down. The OxiRed™ Probe solution is stable for 1 week at 4°C and 1 month at -20°C.

HRP: Dissolve in 220 μl assay buffer, pipetting up and down. The HRP solution is stable for 1 week at 4°C and 1 month at -20°C.

V. Hydrogen Peroxide Assay:

1. Sample Preparations:

Collect cell culture supernatant, serum, plasma, urine and other biological fluids (contains 0.8-6 μM H_2O_2). Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Remove particulate pellet. Assay immediately or aliquot and store the samples at -80°C . Avoid repeated freeze-thaw cycles. Add 2-50 μl samples into each well, bring the volume to 50 μl with assay buffer.

2. H_2O_2 Standard Curve:

For the Colorimetric Assay: Dilute 10 μl 0.88M H_2O_2 standard into 870 μl dH_2O to generate 10 mM H_2O_2 standard, then dilute 10 μl 10 mM H_2O_2 standard into 990 μl dH_2O to generate 0.1 mM H_2O_2 standard. Add 0, 10, 20, 30, 40, 50 μl of the 0.1 mM H_2O_2 standard into 96-well plate in duplicate to generate 0, 1, 2, 3, 4, 5 nmol/well H_2O_2 standard.

For the Fluorometric Assay: Dilute 100 μl of the 0.1 mM H_2O_2 standard into 900 μl dH_2O to generate 10 μM H_2O_2 Standard. Add 0, 10, 20, 30, 40, 50 μl of the 10 μM H_2O_2 standard into 96-well plate in duplicate to generate 0, 0.1, 0.2, 0.3, 0.4, 0.5 nmol/well H_2O_2 standard.

3. Reaction Mix: Mix enough reagents for the number of assays to be performed. For each well, prepare a total 50 μl Reaction Mix:

Colorimetric Assay Fluorometric Assay

46 μl Assay Buffer 48 μl Assay Buffer

2 μl OxiRed™ Probe solution 1 μl OxiRed™ Probe solution

2 μl HRP solution 1 μl HRP solution

Add 50 μl of the Reaction Mix to each test samples and H_2O_2 standards. Mix well.

Incubate at room temperature for 10 min.



*For a more sensitive assay, you can dilute the standard 10 fold further, decrease OxiRed™ amount to 0.2 μl and HRP amount to 0.4 μl per well, it will decrease the fluorescence background and detects as low as 2 pmol/well (or 40 μM concentration) H_2O_2 .

4. Measure OD(570 nm) or fluorescence (Ex/Em = 535/587 nm) in a micro-plate reader.

5. Calculation: Correct background by subtracting the value derived from the 0 nmol H_2O_2 control from all sample and standard readings (Note: The background reading can be significant and must be subtracted from sample readings). Plot the H_2O_2 standard curve. Apply your sample readings to the standard curve. H_2O_2 concentrations of the test samples can then be calculated:

$C = \text{Sa}/\text{Sv}$ (pmol/ μl or μM),

where Sa is the sample amount from your standard curve (in pmol), Sv is sample volume (μl).



VI. Related Products:

Glutathione Reductase Assay Kit
 Glutathione Peroxidase Assay Kit
 Colorimetric Glutathione Detection Kit
 ApoGSH Glutathione Detection Kit
 Glutathione Kit (GSH, GSSG and Total)
 GST Fluorometric Assay Kit
 GST Colorimetric Assay Kit
 Triglyceride Assay Kit
 Acid Phosphatase Assay Kit
 ADP/ATP Ratio Assay Kit
 Phosphate Fluorescence Assay Kit
 Phosphate Colorimetric Assay Kit
 NAD/NADH Quantification Kit
 NADP/NADPH Quantitation Kit
 Pyruvate Assay Kit
 Lactate Assay Kit/ II
 Ammonia Assay Kit
 Glutamate Assay Kit
 Glucose Assay Kit
 Fatty Acid Assay Kit

GENERAL TROUBLESHOOTING GUIDE

Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> • Use of ice-cold assay buffer • Omission of a step in the protocol • Plate read at incorrect wavelength • Use of a different 96-well plate 	<ul style="list-style-type: none"> • Assay buffer must be at room temperature • Refer and follow the data sheet precisely • Check the wavelength in the data sheet and the filter settings of the instrument • Fluorescence: Black plates (clear bottoms) ; Luminescence: White plates ; Colorimeters: Clear plates
Samples with erratic readings	<ul style="list-style-type: none"> • Use of an incompatible sample type • Samples prepared in a different buffer • Cell/ tissue samples were not completely homogenized • Samples used after multiple free-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples 	<ul style="list-style-type: none"> • Refer data sheet for details about incompatible samples • Use the assay buffer provided in the kit or refer data sheet for instructions • Use Dounce homogenizer (increase the number of strokes); observe for lysis under microscope • Aliquot and freeze samples if needed to use multiple times • Troubleshoot if needed • Use fresh samples or store at correct temperatures until use
Lower/ Higher readings in Samples and Standards	<ul style="list-style-type: none"> • Improperly thawed components • Use of expired kit or improperly stored reagents • Allowing the reagents to sit for extended times on ice • Incorrect incubation times or temperatures • Incorrect volumes used 	<ul style="list-style-type: none"> • Thaw all components completely and mix gently before use • Always check the expiry date and store the components appropriately • Always thaw and prepare fresh reaction mix before use • Refer datasheet & verify correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly
Readings do not follow a linear pattern for Standard curve	<ul style="list-style-type: none"> • Use of partially thawed components • Pipetting errors in the reaction mix • Air bubbles formed in well • Calculation errors • Substituting reagents from older kits/ lots 	<ul style="list-style-type: none"> • Thaw and resuspend all components before preparing the reaction mix • Prepare a master reaction mix whenever possible • Pipette gently against the wall of the tubes • Recheck calculations after referring the data sheet • Use fresh components from the same kit
Unanticipated results	<ul style="list-style-type: none"> • Measured at incorrect wavelength • Samples contain interfering substances • Use of incompatible sample type • Sample readings above/below the linear range 	<ul style="list-style-type: none"> • Check the equipment and the filter setting • Troubleshoot if it interferes with the kit • Refer data sheet to check if sample is compatible with the kit or optimization is needed • Concentrate/ Dilute sample so as to be in the linear range
<p>Note# The most probable list of causes is under each problem section. Causes/ Solutions may overlap with other problems.</p>		

GENT