**Dako**  
**EnVision+ System- HRP**  
**Labelled Polymer**  
**Anti-Rabbit**

**Code K4002** 15 mL  
**Code K4003** 110 mL

**ENGLISH**

**Intended use**

For In Vitro diagnostic use.

These instructions apply to the **Dako EnVision®+, Peroxidase** (DC EnVision+ System, HRP).

This kit is intended for use with primary antibodies from **rabbit** supplied by the user for the qualitative identification of antigens by light microscopy in normal and pathological paraffin-embedded tissues, cryostat tissues or cell preparations. Tissues processed in a variety of fixatives including ethanol, B-5, Bouin's, zinc formalin, and neutral buffered formalin may be used.

Refer to the "General Instructions for Immunohistochemical Staining" or the Detection System "Instructions" of IHC procedures for: (1) Principle of Procedure, (2) Materials Required, Not Supplied, (3) Storage, (4) Specimen Preparation, (5) Staining Procedure, (6) Quality Control, (7) Troubleshooting, (8) Interpretation of Staining, (9) General Limitations.

**Summary and explanation**

The EnVision+ System, HRP is a two-step IHC staining technique. This system is based on an HRP labelled polymer which is conjugated with secondary antibodies. The labelled polymer does not contain avidin or biotin. Consequently, nonspecific staining resulting from endogenous avidin-biotin activity in liver, kidney, lymphoid tissues and cryostat sections is eliminated or significantly reduced. The reagent in the EnVision+, HRP is ready-to-use. This system is an extremely sensitive method and, as a result, optimal dilutions of the primary antibody are up to 20 times higher than those used for the traditional PAP technique, and several-fold greater than those used for the traditional ABC or LSAB methods. This protocol offers an enhanced signal generating system for the detection of antigens present in low concentrations or for low titer primary antibodies. Primary antibodies produced in rabbit react well with the labelled polymer. The interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls.

**Principles of procedure**

Any endogenous peroxidase activity is quenched by incubating the specimen with an appropriate endogenous blocking reagent. The specimen is then incubated with an appropriately characterized and diluted rabbit primary antibody, followed by incubation with the labelled polymer, using two sequential 30-minute incubations. ***It should be noted that for antibodies requiring enzyme digestion or target retrieval, it may be necessary to increase incubation times of the primary antibody and labelled polymer by 5 to 10 minutes.*** Staining is completed by a 5-10 minute incubation with a suitable substrate-chromogen.

**Reagent provided**

**Code K4002**

The following is sufficient for 150 tissue sections, based upon 100 µL per section:

<i>Bottle No.</i>	<i>Quantity</i>	<i>Description</i>
2	1x15 mL	<b>Labelled Polymer:</b> Peroxidase labelled polymer conjugated to goat anti-rabbit immunoglobulins in Tris-HCl buffer containing stabilizing protein and an anti-microbial agent.

**Code K4003**

The following is sufficient for 1100 tissue sections, based upon 100 µL per section:

<i>Bottle No.</i>	<i>Quantity</i>	<i>Description</i>
1	1x110 mL	<b>Labelled Polymer:</b> Peroxidase labelled polymer conjugated to goat anti-rabbit immunoglobulins in Tris-HCl buffer containing stabilizing protein and an anti-microbial agent.

**Materials required, but Not Supplied**

Absorbent wipes  
Control tissue, positive and negative  
Counterstain; aqueous based, such as Mayer's Hematoxylin or Lillie's Modified Mayer's Hematoxylin (code S3309)  
Coverslips  
Distilled water  
Endogenous peroxide blocking reagent such as Peroxidase Block (code S2001)  
Ethanol, absolute and 95%  
Light microscope (20x–800x)  
Mounting media, such as GlycerGel® Mounting Medium (code C0563) or Faramount, Aqueous Mounting Medium, Ready-to-use (code S3025) or nonaqueous permanent mounting medium, Ultramount (code S1964)  
Primary antibodies and negative control reagent  
Slides, poly-L-lysine coated or Silanized Slides (code S3003)  
Staining jars or baths  
Substrate-Chromogen solution such as AEC+ Substrate-Chromogen (code K3469, 110 mL, ready-to-use) or Liquid DAB+ (code K3468)  
Timer (capable of 3–40 minute intervals)  
Wash bottles  
Wash Buffer Solution  
Xylene, toluene or xylene substitutes

**Optional materials required, but not supplied**

Ammonium hydroxide, 15 mol/L diluted to 0.037 mol/L  
PAP Pen (code S2002)

**Precautions**

1. For professional users.
2. Do not use reagents beyond expiration date for prescribed storage method. If reagents are stored under any conditions other than those specified in the product insert, they must be validated by the user.
3. Do not substitute reagents from other lot numbers or from kits of other manufacturers.
4. Enzymes and substrate-chromogens may be affected adversely if exposed to excessive light levels. Do not store kit components or perform staining in strong light, such as direct sunlight.
5. Incubation times or temperatures other than those specified may give erroneous results; any such changes must be validated by the user.
6. As with any product derived from biological sources, proper handling procedures should be used.
7. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
8. Unused solution should be disposed of according to local, State and Federal regulations.

**Storage**

EnVision+, HRP is to be stored at 2–8°C. Do not freeze.  
Do not use after expiration printed on reagent vials and product label.

Alteration in the appearance of any reagent, such as precipitation, may indicate instability or deterioration. In such cases, the reagent(s) is (are) not to be used.

There are no obvious signs to indicate instability of these products. Therefore, positive and negative controls should be tested simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variation in laboratory procedures and a problem with the kit is suspected, contact Dako Technical Support.

**Reagent preparation** It is convenient to prepare the following reagents prior to staining.

**Wash Buffer Solution**

TBST, 0.05 mol/L Tris Buffered Saline with Tween (code S3006) is the recommended wash buffer for automated and manual IHC detection. TBS, 0.05 mol/L Tris Buffered Saline (code S1968) and PBS, 0.02 mol/L Phosphate Buffered Saline (code S3024) are also suitable wash buffer solutions for manual staining. Wash buffer solutions containing sodium azide are not recommended. Sodium azide will inactivate horseradish peroxidase (HRP) resulting in negative staining. Store unused buffer at 2–8°C. Discard buffer if cloudy in appearance.

Distilled water may be used for rinsing the peroxidase block, substrate and counterstain.

**Primary Antibody**

NP-Series "Plus" ready-to-use antibodies are optimized and recommended for use with Dako "Plus" high-sensitivity detection systems. Concentrated antibodies are also available from Dako. Optimization of concentrated antibodies is required by the end user. Dilutions should be prepared using Antibody Diluent (code S0809) or a diluent containing 0.05 mol/L Tris-HCl buffer with 1% bovine serum albumin (BSA). Dako N-Series ready-to-use antibodies are not optimized for use with Dako "Plus" detection systems. For most primary antibodies used with this kit, an incubation time of 30 minutes is sufficient.

**Negative Control Reagent**

Ideally, a negative control reagent contains an antibody which exhibits no specific reactivity with human tissues or normal/ nonimmune serum in the same matrix/solution as the diluted primary antibody. The negative control reagent should be the same subclass and animal species as the primary antibody, diluted to the same immunoglobulin or protein concentration as the diluted primary antibody using the same diluent. The incubation period for the negative control reagent should correspond to the primary antibody.

When using Dako NP-Series "Plus" ready-to-use antibodies, Universal Negative Control+ optimized for use with rabbit (code NP001) NP-Series "Plus" ready-to-use antibodies is recommended as a negative control reagent.

For specific information regarding negative control reagent preparation, refer to Quality Control section.

**Counterstain**

The colored end-product of the staining reaction is alcohol insoluble and may be used with alcohol or aqueous-based counterstains such as Mayer's hematoxylin or Lillie's Modified Mayer's Hematoxylin (code S3309). Follow counterstaining of hematoxylin with a thorough rinse in distilled water, then immerse tissue slides into a bath of 0.037 mol/L ammonia or similar bluing agent. 0.037 mol/L ammonia water is prepared by mixing 2.5 mL of 15 mol/L (concentrated) ammonium hydroxide with 1 liter of water.

Unused 0.037 mol/L ammonia may be stored at room temperature (20–25°C) in a tightly capped bottle for up to 12 months.

Consult manufacturers' guidelines for alternative counterstaining procedures.

**Mounting Media**

Glycergel Mounting Medium (code C0563) or Faramount, Aqueous Mounting Medium, Ready-to-use (code S3025) is recommended for aqueous mounting. Glycergel must be heated to at least 50°C just prior to use. Non-aqueous permanent mounting medium (Ultramount, code S1964) may also be used.

**Specimen preparation** Refer to the "General Instructions for Immunohistochemical Staining" and/or the Antibody Specification Sheet.

Prior to IHC staining, tissues must be fixed and processed. Fixation prevents autolysis and putrefaction of excised tissues, preserves antigenicity, enhances the refractive index of tissue constituents and increases the resistance of cellular elements to tissue processing. Tissue processing includes dehydration, clearing of dehydrating agents, infiltration of embedding media, embedding and sectioning of tissues. The most common fixatives for IHC tissue preparations are discussed in the General Instructions for Immunohistochemical Staining. These are guidelines only. Optimal procedures must be determined and verified by the user.

## Staining procedure

### Procedural Notes

The user should read these instructions carefully and become familiar with the product content prior to use. The reagent and instructions supplied have been designed for optimal performance. Further dilution of the reagent or alteration of incubation times or temperatures may give erroneous results. All reagents should be equilibrated to room temperature (20–25°C) prior to immunostaining. Likewise, all incubations should be performed at room temperature.

Do not allow tissue sections to dry during the staining procedure. Dried tissue sections may display increased nonspecific staining. Cover slides exposed to drafts. If prolonged incubations are used, place tissues in a humid environment.

The sensitivity of the EnVision+ System, HRP can be further increased by lengthening the incubation times of Steps 2 and 3 for 5-10 minutes.

### Staining Protocol

#### STEP 1 PEROXIDASE BLOCK

Tap off excess buffer. Using a lintless tissue (such as Kimwipe or gauze pad), carefully wipe around the specimen to remove any remaining liquid and to keep reagent within the prescribed area. Apply enough Peroxidase Block from Bottle 1 to cover specimen. Incubate 5 (±1) minutes. Rinse gently with distilled water or buffer solution from a wash bottle (do not focus flow directly on tissue) and place in a fresh buffer bath.

#### STEP 2 PRIMARY ANTIBODY OR NEGATIVE CONTROL REAGENT

Tap off excess buffer and wipe slides as before. Apply enough optimally diluted primary antibody or negative control reagent to cover specimen. Incubate 30 (±1) minutes. Rinse gently with buffer solution from a wash bottle (do not focus flow directly on tissue) and place in a fresh buffer bath.

If the staining procedure must be interrupted, slides may be kept in a buffer bath following incubation of the primary antibody (Step 2) for up to one hour at room temperature (20–25°C) without affecting the staining performance.

#### STEP 3 PEROXIDASE LABELLED POLYMER

Tap off excess buffer and wipe slides as before. Apply enough Labelled Polymer from Bottle 2 to cover specimen. Incubate 30 (±1) minutes. Rinse slides as in Step 2.

#### STEP 4 SUBSTRATE-CHROMOGEN

Wipe slides as before. Apply enough substrate-chromogen solution to cover specimen. Incubate for 5–10 minutes. Rinse gently with distilled water from a wash bottle (do not focus flow directly on tissue). Collect substrate-chromogen waste in a hazardous materials container for proper disposal.

#### STEP 5 HEMATOXYLIN COUNTERSTAIN (optional)

Immerse slides in a bath of aqueous hematoxylin (code S3309). Length of incubation depends on the strength of hematoxylin used. Rinse gently in a distilled water bath. Dip slides 10 times into a bath of 0.037 mol/L ammonia or similar bluing agent. Rinse slides in a bath of distilled or deionized water for 2–5 minutes.

#### STEP 6 MOUNTING

Specimens may be mounted and coverslipped with an aqueous-based mounting medium such as Glycergel Mounting Medium (code C0563) or Faramount (code S3025) or non-aqueous permanent mounting media, Ultramount (code S1964). Slides can also be dehydrated and permanently mounted.

## Quality control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls.

Consult the quality control guidelines of the College of American Pathologists (CAP) Certification Program for Immunohistochemistry and references 1 through 3 for additional information. Refer to the specification sheet of each primary antibody used for details regarding sensitivity and immunoreactivity.

Refer to the "General Instructions for Immunohistochemical Staining" for further information on positive and negative controls.

## Staining interpretation

Refer to the "General Instructions for Immunohistochemical Staining" for interpretation guidelines.

**Product specific limitations**

Tissue staining is dependent on the proper handling and processing of tissues prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false-negative results.

Use of old or unbuffered fixatives, or exposure of tissues to excessive heat (greater than 60°C) during processing may result in decreased staining sensitivity.

Endogenous peroxidase or pseudoperoxidase activity can be found in hemoproteins such as hemoglobin, myoglobin, cytochrome, and catalase as well as in eosinophils.<sup>4,5</sup> This activity can be inhibited by incubating specimens with Peroxidase Block for five minutes prior to the application of the primary antibody. Blood and bone marrow smears and frozen tissues can also be treated with this reagent. However, this procedure does not abolish the reddish-brown pigment of hemoproteins. Alternately, a solution of methanol-hydrogen peroxide can be used. Some antigens may become denatured with this procedure.

Tissues from persons infected with hepatitis B virus and containing hepatitis B surface antigen (HBsAg) may exhibit nonspecific staining with horseradish peroxidase.<sup>6</sup>

Normal/nonimmune sera from the same animal source as the secondary antisera used in blocking steps may cause false-negative or false-positive results due to auto-antibodies or natural antibodies. The reagents supplied in this kit have been optimally diluted. Further dilution may result in loss of antigen detection.

**Troubleshooting**

Problem	Probable Cause	Suggested Action
1. No staining of any slides.	1a.Reagents not used in proper order.	1a.Review application of reagents.
	1b.Sodium azide in buffer bath.	1b.Use fresh, azide-free buffer.
2. Weak staining of all slides.	2a.Sections retain too much solution after wash bath.	2a.Gently tap off excess solution before wiping around section.
	2b.Slides not incubated long enough with antibodies or substrate mixture.	2b.Review recommended incubation times.
3. Excessive background staining in all slides.	3a.Specimens contain high endogenous peroxidase activity.	3a.Use longer incubation time of peroxidase block, Bottle 1.
	3b.Paraffin incompletely removed.	3b.Use fresh xylene or toluene baths. If several slides are stained simultaneously, the second xylene bath should contain fresh xylene.
	3c.Slides not properly rinsed.	3c.Use fresh solutions in buffer baths and wash bottles. Alternately, use TBST 0.05 mol/L Tris Buffer containing 0.3 mol/L NaCl and 0.1% Tween (code S3306).
	3d.Faster than normal substrate reaction due to e.g. excessive room temperature.	3d.Use shorter incubation time with substrate-chromogen solution.
	3e.Sections dried during staining procedure.	3e.Use humidity chamber. Wipe only three to four slides at a time before applying reagent.
	3f.Nonspecific binding of reagents to tissue section.	3f.Apply a blocking solution containing and irrelevant protein.
	3g.Antibody too concentrated.	3g.see higher dilution of the primary antibody.

NOTE: If the problem cannot be attributed to any of the above causes, or if the suggested corrective action fails to resolve the problem, please call Dako Technical Support for further assistance.

Additional information on staining techniques and specimen preparation can be found in the *Handbook - Immunochemical Staining Methods*<sup>7</sup> (available from Dako), *Atlas of Immunohistology*<sup>8</sup> and *Immunoperoxidase Techniques, A Practical Approach to Tumor Diagnosis*.<sup>9</sup>

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**References**  
**Références**  
**Literatur**

1. Elias JM, et al. Special report: Quality control in immunohistochemistry. *Amer J Clin Pathol* 1989; 92:836
2. National Committee for Clinical Laboratory Standards. Internal quality control testing: principles and definitions; approved guideline. Villanova, PA 1991; Order code C24-A:4
3. Herman GE and Elfont EA. The taming of immunohistochemistry: The new era of quality control. *Biotech & Histochem* 1991; 66:194
4. Escribano LM, et al. Endogenous peroxidase activity in human cutaneous and adenoidal mast cells. *J Histochem Cytochem* 1987; 35:213
5. Elias JM. *Immunohistopathology: A practical approach to diagnosis*. Chicago: Amer Soc of Clin Pathol Press 1990; 46
6. Omata M, et al. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: A possible source of error in immunohistochemistry. *Amer J Pathol* 1980; 73:626
7. Naish SJ (ed). *Handbook-immunochemical staining methods*. Carpinteria: Dako 1989
8. Tubbs RR, et al. *Atlas of immunohistology*. Chicago: Amer Soc Clin Pathol Press 1986
9. Nadji M and Morales AR. *Immunoperoxidase techniques, a practical approach to tumor diagnosis*. Chicago: Amer Soc Clin Pathol Press 1986

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- Cartun RW. Immunohistochemistry of infectious diseases. *J Histotechnol* 1995; 18(3):195
- Heras A, et al. Enhanced labelled-polymer system for immunohistochemistry. XVth Eur Cong Pathol. Copenhagen, Denmark 1995; Sept 3-8
- Bisgaard K and Pluzek K-P. Use of polymer conjugates in immunohistochemistry: A comparative study of a traditional staining method to a staining method utilizing polymer conjugates. Abstract. XXI Intl Cong Intl Acad Pathol and 12th World Cong Acad Environ Pathol. Budapest, Hungary 1996; Oct 20-25
- Pileri SA, et al. EnVision Plus: A new powerful tool for diagnosis and research. Symposium. XXI Intl Cong Intl Acad Pathol and 12th World Cong Acad Environ Pathol. Budapest, Hungary 1996; Oct 20-25
- Bisgaard K. EnVision Plus—Introduction to a new technology. Abstract. Dako Symposium. XXI Intl Cong Intl Acad Pathol. Budapest, Hungary 1996; Oct 20-25
-