**Dako REAL™ Detection System, Alkaline Phosphatase/RED,  
Rabbit/Mouse**

Code K5005

4th edition/ 4ème édition/ 4. Ausgabe

For use with Dako automated immunostaining instruments.  
The kit contains reagents sufficient for 500 tests.

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Page/ Page/ Seite

## ENGLISH

Intended Use .....	4
Summary and Explanation .....	4
Reagents .....	5
A. Materials provided.....	5
B. Materials required but not provided.....	5
Precautions.....	6
Storage .....	7
Reagent Preparation.....	7
Specimen Collection and Preparation .....	7
Procedure .....	8
Quality Control .....	9
Interpretation of Results.....	9
Explanation of symbols .....	22

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### **Intended Use**

For in vitro diagnostic use.

Dako REAL™ Detection System, Alkaline Phosphatase/RED, Rabbit/Mouse, is intended for use in immunohistochemistry together with Dako automated immunostaining instruments. The kit is based on the LSAB (labeled streptavidin-biotin) method and is employed in a three-step procedure. The first step is incubation of the tissue with an optimally diluted primary rabbit or mouse antibody, The second step is incubation with Dako REAL™ Link, Biotinylated Secondary Antibodies (AB2), and the third step is incubation with Dako REAL™ Streptavidin Alkaline Phosphatase (AP). The reaction is visualized by a RED chromogen, which is also included in the kit. To ensure a proper staining it is recommended to use the kit in conjunction with other recommended reagents – listed in the section “Materials required but not provided”.

### **Summary and Explanation**

Dako REAL™ Detection System, Alkaline Phosphatase/RED, Rabbit/Mouse, has been designed to give an optimal staining on Dako automated immunostaining instruments, when using the templates/protocols advised in this package insert.

Prior to staining some routinely fixed, paraffin-embedded tissue sections should be subjected to heat-induced epitope retrieval (HIER) using the Target Retrieval Solution specified in the package insert for the primary antibody. Some primary antibodies require enzymatic pre-treatment of tissue for optimal staining instead of HIER. Please note that Dako REAL™ Proteinase K, code S2019, is designed to be used in combination with HIER in order to ease batch processing of the slides.

Endogenous alkaline phosphatase should be blocked by adding Levamisole (Bottle G) to the AP Substrate Buffer (Bottle F) prior to adding Chromogen Red 1, Red 2 and Red 3. Alternatively, the endogenous alkaline phosphatase can be blocked with Dako Dual Endogenous Enzyme Block, code S2003.

Due to an effective washing procedure and the presence of stabilizing proteins in the Dako reagents, extra blocking steps to reduce non-specific background staining are unnecessary.

Dako Wash Buffer specific for the individual Dako instruments is recommended.

Primary antibodies are not provided with the kit. We recommend the use of Dako N-Series Ready-to-Use Primary Antibodies or Dako concentrated Primary Antibodies.

Dako REAL™ Detection System, Alkaline Phosphatase/RED, Rabbit/Mouse, is based on an indirect streptavidin-biotin method, which has been shown to provide an increased sensitivity compared with analogous methods. The biotinylated secondary antibody reagent of the kit reacts equally well with rabbit and mouse immunoglobulins, thus only one secondary reagent is required for rabbit or mouse primary antibodies. The secondary antibody reagent has been optimally labeled with biotin using a 7-atom spacer arm. This allows each biotinylated antibody molecule to react with several alkaline phosphatase-conjugated streptavidin molecules and thereby increases the sensitivity of the detection system. The alkaline phosphatase-conjugated streptavidin of the kit is prepared from streptavidin isolated from *Streptomyces avidinii*, and alkaline phosphatase isolated from calf intestinal mucosa. The conjugation is performed by a refined two-step glutaraldehyde method. Because the isoelectric point of streptavidin is close to neutral, non-specific ionic interactions with specimen components are avoided. The biotinylated secondary antibody reagent as well as the alkaline phosphatase-conjugated streptavidin is provided as ready-to-use in dropper bottles and have to be poured into appropriate instrument containers before use.

The substrate system consists of four components: Chromogen Red 1, Chromogen Red 2, Chromogen Red 3, and an AP Substrate Buffer. The chromogen is a Fast Red-type chromogen. In addition concentrated Levamisole is provided to be added to the Substrate Working Solution if blocking of endogenous alkaline phosphatase is considered necessary. Before use, the Chromogen Red 1, Red 2 and Red 3 must be diluted in the AP Substrate Buffer. The substrate-chromogen produces a crisp red end product at the site of the target antigen.

Dako REAL™ Hematoxylin, code S2020, or Hematoxylin for the Dako Autostainer, code S3301 is recommended for counterstaining as they provide a clear blue, nuclear staining.

The stained tissue sections may be mounted with either aqueous or organic-solvent-based mounting medium.

## Reagents

### A. Materials provided

Bottle A

LINK, BIOTINYLATED SECONDARY  
ANTIBODIES  
(AB2)

**Dako REAL™ Link, Biotinylated Secondary Antibodies (AB2)**

100 mL, ready-to-use

Biotinylated goat anti-mouse and anti-rabbit immunoglobulins. In buffered solution containing stabilizing protein and sodium azide. The color of this reagent may vary from strong yellow to colorless without having any influence on the performance of the kit.

Bottle B

STREPTAVIDIN  
ALKALINE PHOSPHATASE  
(AP)

**Dako REAL™ Streptavidin Alkaline Phosphatase (AP)**

100 mL, ready-to-use

Streptavidin conjugated to alkaline phosphatase. In buffered solution containing stabilizing protein and preservative.

Bottle C

CHROMOGEN RED 1 (X 28)

**Dako REAL™ Chromogen Red 1**

8 mL, 28x concentrated

Bottle D

CHROMOGEN RED 2 (X 28)

**Dako REAL™ Chromogen Red 2**

8 mL, 28x concentrated

Bottle E

CHROMOGEN RED 3 (X 28)

**Dako REAL™ Chromogen Red 3**

8 mL, 28x concentrated

Bottle F

AP SUBSTRATE BUFFER

**Dako REAL™ AP Substrate Buffer**

250 mL, ready-to-use

Buffered solution containing preservative

Bottle G

LEVAMISOLE (X 501)

**Dako REAL™ Levamisole**

1 mL, 501x concentrated

Buffered solution containing levamisole

### B. Materials required but not provided

Dako automated immunostaining instruments

Dako Wash Buffer specific for the individual instrument

Dual Endogenous Enzyme Block, code S2003 (optional)

Dako REAL™ Proteinase K, code S2019, diluted in Dako REAL™ Proteinase K Diluent, code S2032, or Proteinase K (Ready-to-Use), code S3020 (if necessary)

N-Series Ready-to-Use Primary Antibodies or suitably diluted concentrated primary rabbit or mouse antibodies from Dako

Dako REAL™ Antibody Diluent, code S2022 (for dilution of concentrated antibodies)  
Dako REAL™ Hematoxylin, code S2020 or Hematoxylin for the Dako Autostainer, code S3301  
Target Retrieval Solution. See package insert of the primary antibody  
N-Series Ready-to-Use Rabbit or Mouse Universal Negative control or suitable negative control reagent for the primary antibody  
Dako REAL™ Incubation Container, code S2030 (optional)  
Dako REAL™ Slide Holder, code S2029 (optional)  
Dako instrumentation utensils  
Microscope slides  
General laboratory reagents for dewaxing paraffin-embedded tissue sections  
Mounting medium (aqueous or organic-solvent-based) and coverslips.

## Precautions

1. For professional users.
2. Bottle A, Dako REAL™ Link, Biotinylated Secondary Antibodies (AB2) and Bottle B, Dako REAL™ Streptavidin Alkaline Phosphatase (AP), contain material of mammalian origin. As with any product derived from biological sources, proper handling procedures should be used.
3. Bottle A, Dako REAL™ Link, Biotinylated Secondary Antibodies (AB2), contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. The prepared Chromogen Working Solution must be applied within 20 minutes.
5. Do not expose Bottle C, Dako REAL™ Chromogen Red 1, Bottle D, Dako REAL™ Chromogen Red 2, Bottle E, Dako REAL™ Chromogen Red 3, Bottle G, Dako REAL™ Levamisole, or the Substrate Working Solution (CHROM) to strong light.
6. Bottle E, Dako REAL™ Chromogen Red 3, is very sensitive to heat and exposure to room temperature should therefore be as brief as possible.
7. Bottle E, Dako REAL™ Chromogen Red 3, contains 1–<5% Fast Red KL Salt, and is labeled: Toxic.  
R45 May cause cancer.  
S35 This material and its container must be disposed of in a safe way.  
S45 In case of accident or if you feel unwell, seek medical advice immediately (show label where possible).  
S53 Avoid exposure – obtain special instructions before use.  
Restricted to professional users.  
Caution. This product contains a substance which has not yet been fully tested.  
The classification has been made solely on the basis of a comparison with substances with a similar structure.  
As a main rule, persons under 18 years of age are not allowed to work with this product. Users must be carefully instructed in the proper work procedure, the dangerous properties of the product and the necessary safety instructions. Please refer to the Material Safety Data Sheet (MSDS) for additional information.
8. Bottle G, Dako REAL™ Levamisole, contains 1–<5 % levamisole hydrochloride, and is labeled: Harmful.  
R22 Harmful if swallowed.  
S60 This material and/or its container must be disposed of as hazardous waste.
9. Bottle F, Dako REAL™ AP Substrate Buffer, contains 1–<20 % trometamol. At product concentrations, trometamol does not require hazard labeling. Material Safety Data Sheet (MSDS) available for professional users on request.
10. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
11. Unused solution should be disposed of according to local, State and Federal regulations.

## Storage

Store the kit at 2–8 °C. Do not use after expiration date stamped on the kit. Do not interchange kit components from different lots. If reagents are stored under any conditions other than those specified, the user must verify the conditions. The prepared AP Substrate Buffer with Levamisole should be stored at 2–8 °C in the dark and used within 5 days. The Substrate Working Solution (CHROM) must be used within 20 minutes. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the product is suspected, contact our Technical Services.

## Reagent Preparation

### Substrate Working Solution (CHROM)

The Substrate Working Solution (CHROM) is prepared by mixing thoroughly 750 µL AP Substrate Buffer (Bottle F) with 30 µL Chromogen Red 1 (Bottle C), 30 µL Chromogen Red 2 (Bottle D) and 30 µL Chromogen Red 3 (Bottle E) in that exact order and with thorough mixing after the addition of each chromogen. **Use CHROM within 20 minutes.** This means that the Substrate Working Solution should be prepared and placed on the Dako automated immunostaining instrument no more than 20 minutes before the instrument applies the CHROM on the slides. Discard any volume of CHROM that is not used.

If endogenous alkaline phosphatase is expected to give unwanted unspecific staining, blocking can be achieved by adding Levamisole to the AP Substrate Buffer. Use Levamisole only if required. For those specimens that do require blocking of endogenous alkaline phosphatase, add 1 drop of Levamisole (Bottle G) to each 10 mL of AP Substrate Buffer (Bottle F). Add and mix before adding the Chromogen Red components. Use the AP Substrate Buffer with Levamisole within 5 days (store in the dark at 2–8 °C).

For the preparation of reagents not provided with the kit, please refer to the instructions included with the individual reagent.

## Specimen Collection and Preparation

The specimens may be routinely fixed, paraffin embedded tissue sections and frozen tissue sections.

A variety of fixatives may be used for specimen preservation. Specimens should not be fixed in neutral buffered formalin for more than 24 hours and preferably for a shorter time.

The optimal thickness of paraffin-embedded sections is approximately 4 µm.

The optimal thickness of frozen sections is 4–6 µm.

The specimens should be mounted on microscope slides. The sections should be mounted on the slides as flat and wrinkle-free as possible. Too many wrinkles will have an impact on the staining results.

**NOTE:** The microscope slides must have a width suitable for the Dako instrument. Please refer to the Operator's Manual for the individual Dako instrument for definition of usable microscope slides.

Paraffin sections should be mounted from a pre-heated water bath containing distilled or deionized water. The water bath should contain no additives (such as gelatin, poly-L-lysine etc.). Sections should be dried by heating, generally at a temperature not above 60 °C for a minimum of 60 minutes (e.g. overnight). To ensure proper adherence of sections to slides it is important to drain the water from beneath the sections prior to the oven drying process.

## Procedure

Dako automated immunostaining instruments use techniques based on different principles to obtain an optimal staining result. Before running protocols on the Dako automated immunostaining instrument, please read carefully the **Operator's Manual for the dedicated Dako instrument**.

After deparaffinization and hydration to buffer (water), most tissue sections should need to be subjected to HIER. It is recommended to standardize the HIER in order to ascertain reproducible epitope retrieval. To standardize the treatment, it is recommended to load the slide holder fully with slides, even if only one slide may hold tissue; this ensures an identical heating of sections in every run. At the termination of the epitope retrieval step, slides must be left in the buffer for at least 20 minutes at room temperature. During this cooling, the Dako automated immunostaining instrument can be prepared for the staining run.

A few epitopes do not tolerate the HIER. Please refer to the package insert for the individual Dako primary antibody for instructions.

### *Incubation with primary antibody*

Dako N-Series Ready-to-Use Primary Antibodies can be used with Dako REAL™ Detection System, Alkaline Phosphatase/RED, Rabbit/Mouse, applied on formalin-fixed, paraffin-embedded sections.

A dilution guideline for Dako concentrated primary antibodies and the use of proteinase K pre-treatment in connection with Dako REAL™ Detection System, Alkaline Phosphatase/RED, Rabbit/Mouse, applied on formalin-fixed, paraffin-embedded sections is provided in the package insert of the concentrated primary antibody.

The concentrated primary antibodies should be diluted in Dako REAL™ Antibody Diluent, code S2022.

The slides are loaded onto a Dako instrument and stained using the dedicated template/protocol for Dako REAL™ Detection System, Alkaline Phosphatase/RED, Rabbit/Mouse, code K5005. If the template/protocol is not present on your Dako instrument, please contact our Technical Services.

Before using Dako REAL™ Detection System, Alkaline Phosphatase/RED, Rabbit/Mouse, on the Dako Autostainer for the first time it is necessary to add the reagents and times listed in table 1 to the reagent list of the instrument.

Table 1. Necessary additions to the reagent list and optimal incubation times

Reagent	Code	Suggestion for short name	Template step	Incubation Time
Dual Endogenous Enzyme Block (optional)	S2003	DEEB	End. Enz. Block	5 minutes
Dako REAL™ Proteinase K Solution or Dako Proteinase K	S2019 + S2032 S2030	Pro K	Pre-treatment under antibody	10 minutes
Primary antibodies	Individual	-	Antibody	10*-25 minutes
Dako REAL™ Detection System, Alkaline Phosphatase/RED, Link, Biotinylated Secondary Antibodies (AB2)	K5005	K5005LIN	Secondary Reagent	15 minutes
Dako REAL™ Detection System, Alkaline Phosphatase/RED Strepavidin Alkaline Phosphatase (AP)	K5005	K5005AP	Tertiary Reagent	15 minutes
Dako REAL™ Detection System, Chromogen (RED)	K5005	K5005LPR	Substrate-Batch	20 minutes
Dako REAL™ Hematoxylin Or Dako Hematoxylin for the Dako Autostainer	S2020 S3301	-	Auxiliary	1-2 minutes 1-5 minutes
Blueing Buffer (= Autostainer Buffer)	S 3006	-	Auxiliary	5 minutes

\* N-Series Ready-to-Use Primary Antibodies

For other Dako instruments the reagent names and incubation times are included in the templates/protocols.

## **Quality Control**

Each staining run should include a known positive control specimen to ascertain a proper performance of all the applied reagents. If the positive control specimen fails to demonstrate positive staining, labeling of test specimens should be considered invalid.

A negative control reagent should be used with each specimen to identify any non-specific staining. If non-specific staining cannot be clearly differentiated from the specific staining, the labeling of the test specimen should be considered invalid.

Negative Control for N-Series Rabbit Primary Antibodies, code N1699, or Negative Control for N-Series Mouse Primary Antibodies, code N1698, is recommended as negative control reagent.

## **Interpretation of Results**

The Substrate Working Solution (CHROM) gives a red color at the site of the target antigen recognized by the primary antibody. The red color should be present on the positive control specimen at the expected localization of the target antigen. If non-specific staining is present, this will be recognized as a rather diffuse, red staining on the slides treated with the negative control reagent. Nuclei will be stained blue by the hematoxylin counterstain.

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