

Aconitase Activity Assay Kit

(Catalog #K716-100; 100 reactions; Store kit at 4°C)

I. Introduction:

Aconitase (aconitate hydratase; EC 4.2.1.3) is an iron-sulfur protein containing an $[\text{Fe}_4\text{S}_4]^{2+}$ cluster that catalyzes the stereospecific isomerization of citrate to isocitrate via cis-aconitate in the tricarboxylic acid cycle, a non-redox-active process. Tissue contains two aconitases, a mitochondrial (m-) and a cytosolic (c-) aconitase. They are related, but distinctly different enzymes and are coded for on different chromosomes. Loss of aconitase activity in cells or other biological samples treated with pro-oxidants has been interpreted as a measure of oxidative damage. BioVision's Aconitase Assay Kit is a highly sensitive, simple, direct and HTS-ready colorimetric assay for measuring Aconitase activity in biological samples. In the assay, citrate is converted by aconitase into isocitrate, which is further processed resulting in a product that converts a nearly colorless probe into an intensely colored form ($\lambda = 450\text{nm}$).

II. Kit Contents:

Component	K716-100	Cap Color	Part Number
Assay Buffer	30 ml	WM	K716-100-1
Substrate (lyophilized)	1 vial	Blue	K716-100-2
Developer (lyophilized)	1 vial	Purple	K716-100-3
Enzyme Mix	200 μl	Green	K716-100-4
Cysteine (lyophilized)	1 vial	Red	K716-100-5
$(\text{NH}_4)\text{Fe}(\text{SO}_4)_2$ (lyophilized)	1 vial	Brown	K716-100-6
Isocitrate Standard (10 μmol ; lyophilized)	1 vial	Yellow	K716-100-7

III. Storage and Handling:

Store Enzyme Mix at 4°C and rest of the kit components at -20°C, protected from light. Warm the assay buffer to room temperature before use. Briefly centrifuge vials before opening. Read the entire protocol before performing the assay.

IV. Reagent Reconstitution and General Consideration:

Substrate: Dissolve with 220 μl ddH₂O; sufficient for 100 assays.

Developer: Dissolve with 1.1 ml Assay Buffer before use; sufficient for 100 assays.

Aconitase Activation Solution: Dissolve cysteine and $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2$ with 0.5 ml Assay Buffer separately, and store at -20°C. Take out 0.1 ml cysteine and $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2$ solutions and mix together to prepare fresh activation solution.

All the solutions except the activated Aconitase (store at 4°C) are stable for 1 month at -20°C.

Ensure that the Assay Buffer is at room temperature before use. Keep samples, Enzyme Mix and Aconitase solution on ice during the assay.

V. Aconitase Activity Assay:

1. Sample Preparations:

Homogenize 20 - 40 mg tissue or 10^6 Cells on ice in 0.1 ml cold Assay Buffer; Centrifuge at 800 x g for 10 min at 4°C; Collect the supernatant for c-aconitase assay. For m-aconitase assay, centrifuge the supernatant at 20,000 x g for 15 min at 4°C and collect the pellet, dissolve into 0.1 ml cold Assay Buffer, sonicate for 20 sec. Keep samples at -80°C for storage.

Add 10 μl activation solutions to 100 μl sample; incubate on ice for 1 hr to activate aconitase in the sample.

Add 2 - 50 μl activated samples into each well, and adjust volume to 50 μl . We suggest using a background control group as well as several doses of your sample to ensure the readings are within the linear range.

2. Isocitrate Standard Curve:

Dissolve into 0.5 ml assay buffer to prepare 20 mM isocitrate standard solution. Take 20 μl 20 mM Standard Solution and add 180 μl Assay Buffer to prepare 2 mM isocitrate standard solution. Add 0, 2, 4, 6, 8, 10 μl 2mM Isocitrate Standard solution into 96-well plate in duplicate to generate 0, 4, 8, 12, 16, 20 nmol/well Isocitrate standard. Bring the final volume to 50 μl with Assay Buffer.

3. Reaction Mix: Mix enough reagents for the number of assays to be performed. For each well, prepare a total 50 μl Reaction Mix:

Sample Reaction Mix

46 μl Assay Buffer
2 μl Enzyme Mix
2 μl Substrate

Background Mix

48 μl Assay Buffer
2 μl Enzyme Mix

Add 50 μl of the Sample Reaction Mix to each test samples, background control and Isocitrate standards. Mix well and incubate at 25°C for 30-60 min. Add 10 μl Developer to each well, mix and incubate at 25°C for 10 min. Measure OD 450nm.

4. Calculation: Plot the Isocitrate standard curve. $\Delta\text{OD} = \text{OD}_{\text{sample}} - \text{OD}_{\text{background}}$, apply the ΔOD to the Isocitrate standard curve to get B nmol of isocitrate generated by aconitase in 30 - 60 min.



$$\text{Aconitase Activity} = \frac{B}{T \times V} \times \text{Sample Dilution Factor} = \text{nmol/min/ml} = \text{mU/ml}$$

Where: B is the isocitrate amount from Standard Curve (in nmol)

T is the time incubated (in min)

V is the pretreated sample volume added into the reaction well (in ml)

Unit definition: One unit of Aconitase is the amount of enzyme that will isomerize 1.0 μmol of Citrate to Isocitrate per min at pH 7.4 at 25 °C.



RELATED PRODUCTS:

Colorimetric Glutathione Detection Kit
Glutathione Kit (GSH, GSSG and Total)
GST Colorimetric Assay Kit
Acid Phosphatase Assay Kit
Phosphate Fluorescence Assay Kit
NAD/NADH Quantification Kit
Pyruvate Assay Kit

ApoGSH Glutathione Detection Kit
GST Fluorometric Assay Kit
Triglyceride Assay Kit
ADP/ATP Ratio Assay Kit
Phosphate Colorimetric Assay Kit
NADP/NADPH Quantitation Kit
Lactate Assay Kit/ II

FOR RESEARCH USE ONLY! Not to be used on humans.

GENERAL TROUBLESHOOTING GUIDE:

Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> • Use of ice-cold assay buffer • Omission of a step in the protocol • Plate read at incorrect wavelength • Use of a different 96-well plate 	<ul style="list-style-type: none"> • Assay buffer must be at room temperature • Refer and follow the data sheet precisely • Check the wavelength in the data sheet and the filter settings of the instrument • Fluorescence: Black plates (clear bottoms) ; Luminescence: White plates ; Colorimeters: Clear plates
Samples with erratic readings	<ul style="list-style-type: none"> • Use of an incompatible sample type • Samples prepared in a different buffer • Cell/ tissue samples were not completely homogenized • Samples used after multiple free-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples 	<ul style="list-style-type: none"> • Refer data sheet for details about incompatible samples • Use the assay buffer provided in the kit or refer data sheet for instructions • Use Dounce homogenizer (increase the number of strokes); observe for lysis under microscope • Aliquot and freeze samples if needed to use multiple times • Troubleshoot if needed • Use fresh samples or store at correct temperatures until use
Lower/ Higher readings in Samples and Standards	<ul style="list-style-type: none"> • Improperly thawed components • Use of expired kit or improperly stored reagents • Allowing the reagents to sit for extended times on ice • Incorrect incubation times or temperatures • Incorrect volumes used 	<ul style="list-style-type: none"> • Thaw all components completely and mix gently before use • Always check the expiry date and store the components appropriately • Always thaw and prepare fresh reaction mix before use • Refer datasheet & verify correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly
Readings do not follow a linear pattern for Standard curve	<ul style="list-style-type: none"> • Use of partially thawed components • Pipetting errors in the standard • Pipetting errors in the reaction mix • Air bubbles formed in well • Standard stock is at an incorrect concentration • Calculation errors • Substituting reagents from older kits/ lots 	<ul style="list-style-type: none"> • Thaw and resuspend all components before preparing the reaction mix • Avoid pipetting small volumes • Prepare a master reaction mix whenever possible • Pipette gently against the wall of the tubes • Always refer the dilutions in the data sheet • Recheck calculations after referring the data sheet • Use fresh components from the same kit
Unanticipated results	<ul style="list-style-type: none"> • Measured at incorrect wavelength • Samples contain interfering substances • Use of incompatible sample type • Sample readings above/below the linear range 	<ul style="list-style-type: none"> • Check the equipment and the filter setting • Troubleshoot if it interferes with the kit • Refer data sheet to check if sample is compatible with the kit or optimization is needed • Concentrate/ Dilute sample so as to be in the linear range
<p>Note: The most probable list of causes is under each problem section. Causes/ Solutions may overlap with other problems.</p>		