



DATA SHEET

Technical Data Sheet

Date: 28/4/2006

Lymphosep, Lymphocyte Separation Media CAT N: LM-T1702

Theoretical pH:	7.0 ± 0.5
Colour:	Colourless, clear solution
Storage conditions:	Room temperature
Shelf life:	24 months
Sterility tests:	Bacteria in aerobic and anaerobic conditions Fungi and yeast
Endotoxin:	< 10 EU/ml (<1ng/ml)

Recommended use:

Use aseptic technique when handling this product. For in vitro laboratory use only, not for drug, human or veterinary use.

Application:

Lymphosep is designed for the simple, rapid isolation of lymphocytes from whole blood that has been diluted and treated with anti-coagulant or defibrinating agent.

For best results use blood drawn less than 2 hours before. Do not use blood more than 24 hours from when it was drawn.

Uses:

- 1) Thoroughly mix the Lymphosep by inverting the bottle gently.
- 2) Aseptically transfer 3ml of Lymphosep to a 15ml centrifuge tube
- 3) Mix 2ml of defibrinated or heparinised blood with 2ml of physiological saline (PBS w/o ca w/o Mg) or balanced salt solution (LM-S2041)
- 4) Carefully layer the diluted blood over 3ml of Lymphosep (room temperature) in a 15ml centrifuge, creating a sharp blood-Lymphosep interphase. DO NOT MIX! The quality of the separation is dependant upon a sharp interphase between the lymphocytes and the solution.
- 5) Centrifuge the tube at 400g at room temperature for 15 to 30 minutes. Centrifugation should sediment erythrocytes and polynuclear leukocytes and band mononuclear lymphocytes above the lymphosep.
- 6) Aspirate the top layer of clear plasma to within the 2-3mm above the lymphocyte layer.
- 7) Aspirate the lymphocyte layer plus about half the lymphosep layer below it and transfer it to a centrifuge tube. Add an equal volume of buffered balanced salt solution to the lymphocyte later in the centrifuge tube and centrifuge for 10 minutes at room temperature (18-25°C) at a speed sufficient to sediment to cells without damage i.e. 160-260g. Washing the cells removes Lymphosep and reduces the percentage of platelets.
- 8) Wash the cells again with a buffered balanced salt solution (LM-S2041) and re-suspend in the appropriate medium for your application.

