



## PRODUCT DATA SHEET

**HLA B\*5701 Real-TM**

**Handbook**

**Real Time PCR test for the detection of HLA-B**

**(major histocompatibility complex, class I, B) Allele 5701**

**REF H53-100FRT**

**NAME**

**HLA B\*5701 Real-TM**

**INTRODUCTION**

**Abacavir is a nucleoside reverse-transcriptase inhibitor with activity against the human**

**immunodeficiency virus (HIV), available for once-daily use in combination with other**

**antiretroviral agents, that has shown efficacy, few drug interactions, and a favorable long-term**

**toxicity profile. The most important adverse effect of abacavir that limits its use in therapy and**

**mandates a high degree of clinical vigilance is an immunologically mediated hypersensitivity**

**reaction affecting 5 to 8% of patients during the first 6 weeks of treatment.**

**Symptoms of a**

hypersensitivity reaction to abacavir are nonspecific and include combinations of fever, rash, constitutional symptoms, gastrointestinal tract symptoms, and respiratory symptoms that become more severe with continued dosing. Immediate and permanent discontinuation of abacavir is mandated, resulting in a rapid reversal of symptoms. Subsequent rechallenge with abacavir is contraindicated, since it can result in a more severe, rapid, and potentially life-threatening reaction. In 2002, an association between a diagnosis of hypersensitivity reaction to abacavir and carriage of the major histocompatibility complex class I allele HLA-B\*5701 was reported independently by several independent studies. Studies of cohorts with HIV infection have also shown that avoiding abacavir in HLA-B\*5701 positive patients significantly reduced the incidence of suspected hypersensitivity reaction up to 0,5%. Many clinical studies recommend for this reason, the pharmacogenetic molecular testing of the carriage of the major histocompatibility complex class I allele HLA-B\*5701 in all HIV positive patients treated with abacavir. HLA-B\*5701 Real-TM test can predict who will develop a severe allergic reaction to the anti-

**HIV drug abacavir as the presence of HLA-B\*5701 is significantly associated with an abacavir**

**hypersensitivity.**

#### **INTENDED USE**

**HLA B\*5701 Real-TM is a Real-Time amplification test for the detection of HLA-B (major**

**histocompatibility complex, class I, B) Allele 5701 in the biological materials.**

**The kit HLA B\*5701 Real-TM can be used as screening test for the prevention of abacavir**

**hypersensitivity reactions.**

#### **PRINCIPLE OF ASSAY**

**HLA B\*5701 Real-TM Test is based on two major processes: isolation of genomic DNA from**

**specimens and Real Time amplification with allele specific primers. The real-time PCR**

**monitoring of fluorescence intensities allows the accumulating product detection without**

**reopening of reaction tubes after the PCR run. HLA B\*5701 Real-TM PCR kit is a qualitative**

**test which contains the Internal Control IC (human beta-globine gene), which allows to control**

**the presence of cellular material in the sample.**

#### **MATERIALS PROVIDED**

**Contents Ref. H53-100FRT**

**100 reactions**

**PCR-mix-1 HLA 2 x 0,6 ml**

**PCR-buffer-FRT 2 x 0,3 ml**

**TaqF Polymerase 2 x 0,03 ml**

**Negative Control (C-)\* 4 x 0,5 ml**

**TE-buffer 2 x 0,07 ml**

**Pos C+ (HLA B\*5701 & human DNA ) 1 x 0,2 ml**

**\* must be used in the isolation procedure as Negative Control of Extraction.**

#### **MATERIALS REQUIRED BUT NOT PROVIDED**

**☒ DNA isolation kit**

**☒ Desktop microcentrifuge for ☒eppendorf☒h type tubes**

**☒ Vortex mixer**

**☒ Disposable gloves, powderless**

**☒ Biohazard waste container**

**☒ Refrigerator, Freezer**

**☒ Real Time Thermal cycler**

**☒ Workstation**

**☒ Pipettes (adjustable)**

**☒ Sterile pipette tips with filters**

**☒ Tube racks**

#### **WARNINGS AND PRECAUTIONS**

**1. Wear disposable gloves, laboratory coats and eye protection when handling specimens and**

**reagents. Thoroughly wash hands afterward.**

**2. Do not pipette by mouth.**

- 3. Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.**
- 4. Do not use a kit after its expiration date.**
- 5. Do not mix reagents from different kits.**
- 6. Dispose of all specimens and unused reagents in accordance with local regulations.**
- 7. Specimens and controls should be prepared in a laminar flow hood.**
- 8. Heparin has been shown to inhibit reaction. Use of heparinized specimens is not recommended.**
- 9. Avoid repeated thawing and freezing of the reagents, this may reduce the sensitivity of the test.**
- 10. Once the reagents have been thawed, vortex and centrifuge briefly the tubes.**
- 11. Prepare quickly the Reaction mix on ice or in the cooling block.**
- 12. Specimens may be infectious. Use Universal Precautions when performing the assay.**
- 13. Clean and disinfect all spills of specimens or reagents using a disinfectant such as 0,5% sodium hypochlorite, or other suitable disinfectant. Follow by wiping down the surface with 70% ethanol.**
- 14. Avoid contact of specimens and reagents with the skin, eyes and mucous membranes. If**

these solutions come into contact, rinse immediately with water and seek medical advice

immediately.

15. Material Safety Data Sheets (MSDS) are available on request.

16. Use of this product should be limited to personnel trained in the techniques of DNA

amplification.

17. Workflow in the laboratory must proceed in a uni-directional manner, beginning in the

Extraction Area and moving to the Amplification Area. Do not return samples, equipment

and reagents in the area where you performed previous step. Personnel should be using

proper anti-contamination safeguards when moving between areas.

#### STORAGE INSTRUCTIONS

The kit HLA B\*5701 Real-TM must be stored at -20°C.

The kit can be shipped at 2-8°C for 3-4 days but should be stored at -20°C immediately on

receipt.

#### STABILITY

HLA B\*5701 Real-TM is stable up to the expiration date indicated on the kit label. The product

will maintain performance through the control date printed on the label.

Exposure to light, heat

or humidity may affect the shelf life of some of the kit components and should be avoided.

Repeated thawing and freezing of these reagents should be avoided, as this may reduce the sensitivity.

## **SAMPLE COLLECTION, STORAGE AND TRANSPORT**

HLA B\*5701 Real-TM can analyze genomic DNA extracted from:

☑E whole blood collected in EDTA tubes;

Specimens can be stored at +2-8°C for no longer than 48 hours, or freeze at -20°C to -80°C.

Transportation of clinical specimens must comply with country, federal, state and local

regulations for the transport of etiologic agents.

## **DNA ISOLATION**

The following isolation kits are recommended:

☑E Genomic column DNA Express . spin column extraction kit (Sacace, REF K-1-1/E)

☑E QIAamp DNA Blood mini kit (Qiagen, 51104)

DNA samples with concentration in range from 20 to 200 ng/μl could be analyzed.

## **REAGENT PREPARATION**

1. Prepare the required quantity of reaction tubes for samples (N) and controls (N+2).

2. Prepare Reaction Mix by adding for each sample into the new sterile tube 10\*N μl of PCRmix-

1 HLA and 5\*N μl of PCR-buffer-FRT and 0,5\*N μl of TaqF Polymerase. Vortex and centrifuge for 2-3 sec.

**3. Add to each reaction tube 15 µl of Reaction Mix and 10 µl of extracted genomic DNA. Mix**

**by pipetting.**

**4. Prepare for each panel 2 controls:**

**• add 10 µl of TE-buffer to the tube labeled Amplification Negative Control;**

**• add 10 µl of Pos C+ to the tube labeled Amplification Positive Control;**

**5. Insert the tubes in the thermalcycler.**

**Create a temperature profile on your Real-time instrument as follows:**

**Rotor type instruments1 Plate type or modular instruments2**

**Stage „Semp, 2“, R Time**

**Fluorescence**

**detection**

**Cycle**

**repeats**

**„Semp, 2“, R Time Fluorescence detection**

**Cycle**

**repeats**

**Hold 95 15 min . 1 95 15 min . 1**

**Cycling**

**95 5 s .**

**5**

**95 5 s .**

**5**



**60 20 s . 60 20 s .**

**Cycling 2**

**95 5 s .**

**40**

**95 5 s .**

**40**

**60 40 s**

**FAM(Green),**

**JOE(Yellow)**

**60 50 s FAM, JOE/HEX/Cy3**

**1 For example Rotor-Gene. 3000/6000 (Corbett Research, Australia)**

**2 For example, iQ5./iQ iCycler. (BioRad, USA); Mx3000P/Mx3005P.  
(Stratagene, USA), Applied**

**BiosystemsR 7300/7500 Real Time PCR (Applied), SmartCyclerR (Cepheid)**

**.**

## **RESULTS ANALYSIS**

**The results are interpreted by the device software through the presence of crossing of**

**fluorescence curve with the threshold line.**

**DNA HLA\*B5701 is detected on the JOE (Yellow)/HEX/Cy3 channel and IC on the FAM**

**(Green) channel.**

**Results are accepted as relevant if both positive and negative controls of amplification along**

with negative control of extraction are passed (see table 1).

#### Table 1. Results for controls

Control Stage for control Ct channel Fam (Green) Ct channel Joe

(Yellow)/ HEX/Cy3 Interpretation

NCS DNA isolation Neg Neg Valid result

TE-buffer Amplification Neg Neg Valid result

Pos C+ Amplification Pos Pos Valid result

☒E The sample is considered to be positive if in the channel Joe (Yellow)/HEX/Cy3 the

result is positive and the value of Ct on this channel is higher than Ct on the Fam (Green)

channel but not more than 5 cycles.

☒E The sample is considered to be negative if in the channel Joe (Yellow)/HEX/Cy3 value is

negative or if the value of Ct on this channel is higher than Ct on the Fam (Green) of

more than 5 cycles.

☒E Normal difference between Joe (Yellow) and Fam (Green) Ct values is 2-3 cycles.

#### TROUBLESHOOTING

1. Absent signal of the IC (Fam (Green) channel): retesting of the sample is required.

☒E The PCR was inhibited.

☒E Make sure that you use a recommended DNA extraction method and follow

the manufacturer's instructions.

☒E The reagents storage conditions didn't comply with the instructions.

☒E Check the storage conditions

☒E The PCR conditions didn't comply with the instructions.

☒E Check the PCR conditions and for the IC detection select the fluorescence channel reported in the protocol.

☒E No correct sample collection or preparation.

**2. No signal on the Joe (Yellow)/Cy3/HEX and Fam (Green) channels with Positive Control.**

☒E The reagents storage conditions didn't comply with the instructions.

☒E Check the storage conditions

☒E The PCR conditions didn't comply with the instructions.

☒E Check the temperature profile and select the fluorescence channel reported

in the protocol.

☒E Incorrect configuration of the PCR reaction:

☒E Check the reagents preparation step.

**3. Any signal with Negative Control.**

☒E Contamination during PCR preparation procedure. All samples results are invalid.

☒E Decontaminate all surfaces and instruments with sodium hypochlorite and ethanol or special DNA decontamination reagents.

☒E Pipette the Positive controls at the end.

☒E Repeat the PCR preparation with the new set of reagents.

**4. Variation of more than 5 cycles between Ct values of Fam(Green) and Joe (Yellow) in a sample.**

**☒E Contamination with genomic material during PCR preparation procedure. Sample**

**result is invalid, retesting of sample is required.**

#### **PERFORMANCE CHARACTERISTICS**

**The analytical sensitivity**

**The analytical sensitivity was determined using Company standard sample POS HLA B\*5701**

**which is a plasmid pGEM-t solution containing DNA fragment of allele HLA B\*5701 (2nd and**

**3rd exon of HLA B\*5701 gene) with a known concentration. The sample was tested using real**

**time PCR with its dilutions 10x and 2x in standard DNA buffer.**

**The analytical sensitivity of Sacace HLA B\*5701 Real-TM kit was declared as 1000 copies/ml**

**(positive result in a 12 times tested sample with this concentration).**

**The analytical specificity**

**The analytical specificity was valued using DNA controls (13 samples) with a known allelic**

**combination of HLA previously tested using sequence specific primers.**

#### **DNA CONTROLS**

**No. Sample Allele HLA-B**

**1 R98-149776Z 5704 5301**

**2 R97-403860J 5701 801**

**3 R98-117211A 5701 3701**

**4 R98-106525C 4201 5801**

**5 R98-136122N 5801 1502**

**6 R98-1036872Z 4403 4001**

**7 R98-278853J 15 8**

**8 R98-101966H 5701 801**

**9 R97-397992J 44 40**

**10 R90 50043X 35 61**

**11 R91 7077Y 51 55**

**12 BSN 9400419 5101 5703**

**13 R92 20601H 27 40**

**DNA samples performed on RotorGene 6000**

**Samples containing**

**HLA B\*5701allele**

**Amplification of**

**s-globine gene (IC)**

**The results of the analysis of DNA samples show the absence of a cross reaction with above**

**mentioned alleles, especially belonging to B\*57 group (B\*5703 and B\*5704).**

**REPRODUCIBILITY**

**Reproducibility was tested in two independent assays using different dilutions of company**

**standard sample POS HLA B\*5701**

**Test**

**Concentration,**

**copies/ml**

**Pos/total Medium Ct Dev Ct**

**1**

**2 \* 107 4/4 15,73 0,12**

**2 \* 106 4/4 18,91 0,14**

**2 \* 105 4/4 22,28 0,12**

**2 \* 104 4/4 25,90 0,24**

**2 \* 103 4/4 28,01 0,29**

**1 \* 103 4/4 39,25 0,36**

**2**

**2 \* 107 3/3 15,40 0,11**

**2 \* 106 3/3 18,58 0,19**

**2 \* 105 3/3 22,93 0,18**

**2 \* 104 3/3 25,56 0,24**

**2 \* 103 3/3 28,78 0,27**

**1 \* 103 9/9 29,93 0,37**

**No. Sample Allele HLA-B**

**Results with**

**Sacace. HLA B\*5701 Real-TM kit**

**HLA B\*5701 „B,,K,,O**

**1 R98-149776Z 5704 5301 - +**

- 2 R97-403860J 5701 801 + +  
3 R98-117211A 5701 3701 + +  
4 R98-106525C 4201 5801 - +  
5 R98-136122N 5801 1502 - +  
6 R98-1036872Z 4403 4001 - +  
7 R98-278853J 15 8 - +  
8 R98-101966H 5701 801 + +  
9 R97-397992J 44 40 - +  
10 R90 50043X 35 61 - +  
11 R91 7077Y 51 55 - +  
12 BSN 9400419 5101 5703 - +  
13 R92 20601H 27 40 - +

#### **Diagnostic sensitivity and specificity**

The method of direct sequencing was used as reference because no other commercial certified kit was found.

#### **Results with**

**Sacace. HLA B\*5701 Real-TM kit**

+ -

**Direct**

**sequencing**

**+ 34 0**

**- 0 317**

The experimental group included 351 patients were tested both by direct sequencing method and

by real time PCR using Sacace. HLA B\*5701 Real-TM kit. A complete correspondence was

observed between both methods showing a diagnostic sensitivity and specificity of Sacace HLA

B\*5701 Real-TM kit to 100%.

#### References

☒E High sensitivity of human leukocyte antigen-B\*5701 as a marker of immunologically confirmed

abacavir hypersensitivity in white and black patients. Saag M et al. Clin Infect Dis 46: 1111 .

1118, 2008.

☒E Association between presence of HLA-B\*5701, HLA-DR7, and HLA-DQ3 and hypersensitivity to

HIV-1 reverse-transcriptase inhibitor abacavir. Mallal S, Nolan D et al. Lancet. 2002 Mar

2;359(9308):722-3.

☒E HLA-B\*5701 screening for hypersensitivity to abacavir. Mallal et al N Engl J Med. 2008 Feb

7;358(6):568-79.

☒E Value of the HLA-B\*5701 allele to predict abacavir hypersensitivity in Spaniards. Rodriguez-

Novoa S, Garcia-Gasco P, Blanco F, Gonzalez-Pardo G, Castellares C, Moreno V, Jimenez-

Nacher I, Gonzalez-Lahoz J, Soriano V. AIDS Res Hum Retroviruses. 2007 Nov;23(11):1374-6.



**☒E Prospective HLA-B\*5701 screening and abacavir hypersensitivity: a single centre**

**experience. Waters LJ, Mandalia S, Gazzard B, Nelson M. AIDS. 2007 Nov 30;21(18):2533-4.**

**☒E Abacavir hypersensitivity reaction in primary HIV infection. Stekler J, Maenza J, Stevens C, Holte**

**S, Malhotra U, McElrath MJ, Corey L, Collier AC. AIDS. 2006 Jun 12;20(9):1269-74.**

**☒E The pharmacogenetics of antiretroviral therapy. Phillips EJ. Curr Opin HIV AIDS. 2006**

**May;1(3):249-56.**

#### **EXPLANATION OF SYMBOLS**

**REF Catalogue Number**

**RUO Research Use Only**

**LOT Lot Number**

**Expiration Date**

**Contains reagents**

**Caution!**

**VER Version**

**Manufacturer**

**Temperature limitation**

**\*\*iCycler. and iQ5. are trademarks of Bio-Rad Laboratories**

**\* Rotor-Gene. Technology is a registered trademark of Corbett Research**

**\*MX3000PR and MX3005PR are trademarks of Stratagene**

**\*Applied BiosystemsR is trademarks of Applied Biosystems Corporation**

\* SmartCyclerR is a registered trademark of Cepheid

GENTAUR

Gentaur Molecular Products  
Voortstraat 49  
1910 Kampenhout, Belgium