



---

**Version 2.0808**

**EpiQuik™ *In Situ* Histone H3-K27  
Tri-Methylation Assay Kit**

Catalog No. P-3014T

**User Guide\***

**\*Always use the most updated User  
Guide included in your current order.**

Epigenetek Group Inc.  
760 Parkside Avenue  
Brooklyn, NY 11226  
Tel: 1 877-374-4368  
Fax: 1 718-484-3956  
info@epigenetek.com

---

**www . EPIGENETEK . com**

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

# TABLE OF CONTENTS

---

Introduction . . . . .	3
Principle and Procedure . . . . .	4
Product Use Information . . . . .	5
Kit Contents . . . . .	6
Shipping and Storage . . . . .	6
Materials Required But Not Supplied . . . . .	6
Protocol . . . . .	7
Troubleshooting . . . . .	9
Ordering Information . . . . .	10

Epigentek Group Inc.  
760 Parkside Avenue  
Brooklyn, NY 11226  
Tel: 1 877-374-4368  
Fax: 1 718-484-3956  
info@epigentek.com


## INTRODUCTION

Epigenetic activation or inactivation of genes plays a critical role in many important human diseases, especially in cancer. A major mechanism for epigenetic inactivation of the genes is methylation of CpG islands in genome DNA caused by DNA methyltransferases. Histone methyltransferases (HMTs) control or regulate DNA methylation through chromatin-dependent transcription repression or activation. HMTs transfer 1-3 methyl groups from S-adenosyl-L-methionine to the lysine and arginine residues of histone proteins. G9a and polycomb group enzyme such as EZH2 are histone methyltransferases that catalyze methylation of histone H3 at lysine 27 (H3-K27) in mammalian cells. Tri-methylation of H3-K27 is a facultative heterochromatin mark, which promotes the recruitment of polycomb group proteins for gene silencing. Increased global H3-K27 tri-methylation is also found to be involved in some pathological processes such as cancer progression. There are only a couple of methods such as western blot, which are used for measuring histone H3-K27 tri-methylation. The methods available so far are time consuming and labor intensive, or have low throughput. The *EpiQuik™ In Situ* H3-K27 Tri-Methylation Assay Kit addresses these problems by using a unique procedure to measure *in situ* tri-methylation of histone H3-K27. The kit has the following features:

- Quick and efficient procedure, which can be finished within 3 hours.
- Innovative colorimetric assay without the need for radioactivity, electrophoresis, or chromatography.
- Measurement of *in situ* histone H3-K27 tri-methylation without the need to prepare cell lysates.
- Microplate format makes the assay suitable for high throughput analysis of agents that increases or inhibits H3-K27 tri-methylation.
- Simple, reliable, and consistent assay conditions.

## PRINCIPLE AND PROCEDURE

*EpiQuik™ In Situ* Histone H3-K27 Tri-Methylation Assay Kit is a whole cell-based detection of tri-methylated H3-K27. In this assay, adherent cells are cultured in conventional 96-well microplates. After your experimental treatment, cells are fixed and permeabilized. The tri-methylated H3-K27 is then detected by an anti-tri-methyl H3-K27 antibody. The ratio or amount of tri-methylated H3-K27 can be quantified through HRP conjugated secondary antibody-color development system and is proportional to the intensity of color development.



Schematic Procedure for Using the *EpiQuik™ In Situ* Histone H3-K27 Tri-Methylation Assay Kit

## PRODUCT USE INFORMATION

The *EpiQuik™ In Situ* Histone H3-K27 Tri-Methylation Assay Kit is suitable for specifically measuring histone H3-K27 tri-methylation *in situ* using cultured adherent cells.

The *EpiQuik™ In Situ* Histone H3-K27 Tri-Methylation Assay Kit is for research use only and is not intended for diagnostic or therapeutic application.

Epigentek guarantees the performance of all products in the manner described in our product instructions.

Epigentek reserves the right to change or modify any product at any time to enhance its performance and design.

*EpiQuik™* is a trademark of Epigentek Group Inc.

## KIT CONTENTS

Components	96 assays	2 x 96 assays
	P-3014T-096	P-3014T-192
<b>GK1</b> (10X Wash Buffer)	30 ml	2 x 30 ml
<b>GK2</b> (Permeabilizing Buffer)	30 ml	2 x 30 ml
<b>GK3</b> (Blocking Buffer)	20 ml	2 x 20 ml
<b>GK4</b> (Antibody Buffer)	10 ml	20 ml
<b>GK5</b> (Capture Antibody, 1000 µg/ml)*	9 µl	18 µl
<b>GK6</b> (Detection Antibody, 200 µg/ml)*	10 µl	20 µl
<b>GK7</b> (Developing Solution)	12 ml	24 ml
<b>GK8</b> (Stop Solution)	6 ml	12 ml
30% H <sub>2</sub> O <sub>2</sub> Solution	0.5 ml	1 ml
Tri-Methyl H3-K27 Control (20 µg/ml)	15 µl	30 µl
8-Well Control Strips	2	4
Microplate	1	2
User Guide	1	1

\* For maximum recovery of the products, centrifuge the original vial prior to opening the cap.

## SHIPPING AND STORAGE

The kit is shipped in two parts: one part at ambient room temperature and the second part on frozen ice packs at 4°C.

Upon receipt: (1) Store **GK6** and **Tri-Methyl H3-K27 Control** at -20°C; (2) Store **all other components** at 4°C away from light. The kit is stable for up to 6 months from the shipment date, when stored properly.

**Note:** Check if buffers **GK1** and **GK4** contain salt precipitates before using. If so, warm (at room temperature or 37°C) and shake the buffers until the salts are re-dissolved.

## Materials Required But Not Supplied

Pipettes and pipette tips  
Microplate reader  
15 ml conical tube  
1.5 ml microcentrifuge tubes  
37% formaldehyde  
PBS

## PROTOCOL

Before starting, perform the following:

1. Prepare the following required solutions (not included): **37% formaldehyde**.
2. Ensure that all buffers are in clear solution. Shake or vortex if these buffers precipitate.
  1. Inoculate and grow adherent cells in 96-well microplate to 50-70% confluency. Leave 2-4 wells with no cell inoculation as the blank. Treat cells with appropriate amount of agents that may increase or reduce H3-K27 tri-methylation for appropriate time.
  2. Prepare fixing solution by adding 2.16 ml of **37% formaldehyde** to 18 ml of PBS. Remove culture media from the wells with a wrist-flick.
  3. Immediately add 150  $\mu$ l of fixing solution slowly to the wells and incubate at room temperature for 15 minutes. Remove fixing solution from wells with a wrist-flick; while still inverted, tap the plate gently onto the absorbent paper to remove any excess fixing agent still within the wells.
  4. Dilute **GK1** (10X Wash Buffer) with distilled water (pH 7.2 to 7.5) at a 1:10 ratio (e.g., 1 ml of **GK1** + 9 ml of distilled water). Wash wells once (for 2 minutes) with 150  $\mu$ l of the **diluted GK1**.
  5. Remove wash buffer with a wrist flick; while still inverted, tap the plate onto the absorbent paper to remove any excess solution. Add 150  $\mu$ l of **GK2** to each well and incubate at room temperature for 5 minutes. Meanwhile, prepare 1%  $H_2O_2$  by adding 330  $\mu$ l of **30%  $H_2O_2$  Solution** into 10 ml of **1X GK2**.
  6. Remove **GK2** from the wells with a wrist flick. Add 100  $\mu$ l of the 1%  $H_2O_2$  solution into each well and incubate at room temperature for 10 minutes to remove endogenous peroxidase.
  7. Remove the  $H_2O_2$  solution from wells with a wrist flick and wash the wells twice with 150  $\mu$ l of **diluted GK1**.
  8. Remove wash buffer with a wrist flick; while still inverted, tap the plate onto the absorbent paper to remove any excess solution. Add 150  $\mu$ l of **GK3** to the wells and incubate at 37°C for 45 minutes. Meanwhile, add 50  $\mu$ l of **diluted GK1** to the desired number of control strip wells, followed by adding 1  $\mu$ l of **Tri-Methylated H3-K27 Control** protein at the different amounts (ex: 0.5-20 ng, diluted with distilled water) and incubate at room temperature for 30-45 minutes. For the blank wells, do not add any methylated H3-K4 control protein.
  9. Remove **GK3** with a wrist flick; while still inverted, tap the plate onto the absorbent paper. Wash the wells twice with 150  $\mu$ l of **diluted GK1**. For each wash, remove the wash buffer with a wrist flick; while still inverted, tap the plate onto the absorbent paper to remove any excess solution. Meanwhile, aspirate the solution from the control strip wells and wash the wells with 150  $\mu$ l of **diluted GK1** three times.

10. Dilute **GK5** (at a 1:1000 ratio) to 1  $\mu\text{g/ml}$  with **GK4**. Add 50  $\mu\text{l}$  of **diluted GK5** to the sample wells and *Tri-Methyl H3-K27 control strip wells*. Incubate at room temperature for 60 minutes on an orbital shaker (100 rpm).
11. Remove solution from the wells with a wrist flick and wash the wells four times with 150  $\mu\text{l}$  of **diluted GK1**. For each wash, remove the wash buffer with a wrist flick; while still inverted, tap the plate onto absorbent paper.
12. Dilute **GK6** (at a 1:1000 ratio) to 0.2  $\mu\text{g/ml}$  with **1X GK4**. Add 50  $\mu\text{l}$  of **diluted GK6** to the wells and incubate at room temperature for 30 minutes.
13. Remove solution from the wells with a wrist flick and wash the wells four times with 150  $\mu\text{l}$  of **diluted GK1**. For each wash, remove the wash buffer with a wrist flick; while still inverted, tap the plate onto absorbent paper.
14. Add 100  $\mu\text{l}$  of **GK7** to the wells and incubate at room temperature for 2-10 minutes away from light. Monitor the color development in the sample and control wells (blue).
15. Add 50  $\mu\text{l}$  of **GK8** to the wells and read absorbance on a microplate reader at 450 nm.
16. Calculate % H3-K27 tri-methylation:

$$\text{Methylation \%} = \frac{\text{O.D. (treated sample - blank)}}{\text{O.D. (untreated control - blank)}} \times 100\%$$

17. For accurate calculation, plot OD value versus amount of **Tri-Methylated H3-K27 Control** and determine the slope as  $\Delta \text{OD}/\text{ng}$ .

Calculate tri-methylated H3-K27 amount using the following formula:

$$\text{Tri-Methyl H3-K27 (ng)} = \frac{\text{OD (sample - blank)}}{\text{slope}}$$



## TROUBLESHOOTING

### **No Signal for Both the Positive Control and the Samples**

- |   |  |
|---|--|
| Reagents are added incorrectly.               | Check if reagents are added in order and if any steps of the procedure may have been omitted by mistake. |
| Incubation time and temperature is incorrect. | Ensure the incubation time and temperature described in the protocol are followed correctly.             |

### **No Signal for Only the Sample**

- |   |   |
|---|---|
| Cells are not fixed and permeabilized sufficiently.   | Ensure fixation solution and permeabilizing solution are sufficiently added into cells and incubation time is enough. |
| The protein amount is added into well insufficiently. | Ensure extract contains a sufficient amount of proteins.  |

### **High Background Present for the Blank**

- |                                |  |
|--------------------------------|--|
| The well is not washed enough. | Check if wash at each step is performed according to the protocol. |
| Overdevelopment.               | Decrease development time in step 14.                              |

## ORDERING INFORMATION

<b>Products</b>	<b>Size</b>	<b>Cat. No.</b>
<i>EpiQuik™ In Situ</i> Histone H3-K27 Tri-Methylation Assay Kit	96 assays 2x96 assays	P-3014T-096 P-3014T-192

<b>Available Related Products</b>	<b>Cat. No.</b>
-----------------------------------	-----------------

<i>EpiQuik™</i> DNA Methyltransferase Activity/Inhibition Assay Kit	P-3001
<i>EpiQuik™</i> Histone Methyltransferase Activity/Inhibition Assay Kit (H3-K4)	P-3002
<i>EpiQuik™</i> Histone Methyltransferase Activity/Inhibition Assay Kit (H3-K9)	P-3003
<i>EpiQuik™ In Situ</i> Histone H3-K9 Methylation Assay Kit	P-3016
<i>EpiQuik™</i> Global Histone H3-K4 Methylation Assay Kit	P-3017
<i>EpiQuik™</i> Global Histone H3-K9 Methylation Assay Kit	P-3018
<i>EpiQuik™</i> DNA Demethylase Activity/Inhibition Assay Kit	P-3019
<i>EpiQuik™</i> Global Histone H3-K27 Methylation Assay Kit	P-3020

***Need more components? You can also order parts separately by calling 1-877-374-4368 or e-mailing [sales@epigentek.com](mailto:sales@epigentek.com).***



FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

