**Dako REAL™
Peroxidase-Blocking Solution
Code S2023**

ENGLISH

Intended use	For in vitro diagnostic use. Dako REAL™ Peroxidase-Blocking Solution is intended for use in immunohistochemistry. It is a stable, ready-to-use reagent specially designed for optimal performance when used on tissue sections together with Dako REAL™ Detection Systems on Dako automated immunostaining instruments.
Reagent provided	Phosphate buffer containing hydrogen peroxide, 15 mmol/L NaN ₃ and detergent. Each bottle contains 250 mL ready-to-use reagent.
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃) a chemical highly toxic in pure form. At product concentrations though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. Safety data sheet available for professional users on request. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2–8 °C. Do not use after expiration date stamped on bottle. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services.
Staining procedure	Dako REAL™ Peroxidase-Blocking Solution, Code S2023, is ready-to-use. According to choice any of the Dako REAL™ Detection Systems on the enclosed list may be employed for visualization. These kits provide procedures and package inserts for the primary antibody for heat-induced epitope retrieval on formalin-fixed, paraffin-embedded tissues and procedures for immunostaining using Dako automated immunostaining instruments.
Additional Information	Please refer to package inserts for Dako REAL™ Detection Systems and the individual Operator's Manual of the Dako automated immunostaining instrument for further information on the instrument and instructions for use.

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


Companion products/ Produits Annexes/ Begleitprodukte

Dako REAL™ Detection Systems/ Trousses de détection Dako REAL™/ Dako REAL™ Detektionssysteme
 Alkaline Phosphatase/RED, Rabbit/Mouse, Code K5005
 APAAP, Mouse, Code K5000
 EnVision™ Detection System, Peroxidase/DAB+, Rabbit/Mouse, Code K5007
 Peroxidase/AEC, Rabbit/Mouse, Code K5003
 Peroxidase/DAB+, Rabbit/Mouse, Code K5001

Dako REAL™ Accessory Reagents/ Réactifs auxiliaires Dako REAL™/ Dako REAL™ Zusatzreagenzien
 Antibody Diluent, Code S2022
 Hematoxylin, Code S2020
 Peroxidase-Blocking Solution, Code S2023
 Proteinase K (40x), Code S2019
 Proteinase K Diluent, Code S2032
 Target Retrieval Solution (10x), Code S2031

Dako REAL™ Accessory Utensils/ Ustensiles auxiliaires Dako REAL™/ Dako REAL™ Zubehör
 Incubation Container, Code S2030
 Slide Holder, Code S2029

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