



## **GLYCEROL MOUNTING Medium**

**Catalog number:** AR-6508-02 (100 ml), AR-6508-04 (500 ml), AR-6508-05 (1000 ml).

**Description:** This mounting medium is made with glycerol. It is excellent for preserving thick sections that can not be mounted with mounting medium or fresh frozen tissues that contain a lot of lipids e.g. brain tissue. It is also useful for preserving small embryos. Not recommended for phycoerythrin (PE), phycocyanin (PC) and allophycocyanin (APC)

**Intended Use:** Mounting tissues and cell etc.

**Reagent:** Ready to use mounting medium

**Storage:** Store 2-8°C is recommended, protect from light,

**Refractive Index:** 1.47

**Procedure: Dehydration of specimen is not required.**

1. Bring the vial to room temperature. Turn the vial upside down and open the dropper to remove any air bubbles.
2. Apply 1-3 drops of mounting medium directly on top of the specimen.
3. Place coverslip, it may be necessary to seal the edges with nail polish or any organic mounting medium like Organo mounting medium (IBSC AR-6504)

**Limitation and warranty:** Our warranty is limited to the actual price paid for the product. We are not liable for any property damage, personnel injury, time, effort or economic loss due to our product.

**MSDS:** Avoid skin and eye contact with all laboratory products. Use appropriate laboratory gear, lab coat, gloves and safety glasses. Do not ingest any laboratory products. This product is not approved for administration in human or animals.

---

**For research use only; not for use in diagnostic procedures. FOR IN VITRO LABORATORY USE ONLY**  
“*In vitro* laboratory products for research”



Gentaaur Molecular Products  
Voortstraat 49  
1910 Kampenhout, Belgium