

BM-2 Terrestrial Orchid Medium w/out Agar

Properties

Form:	Powder
Appearance:	White to Cream
Application:	Plant Tissue Culture
Solubility:	Soluble in Water
Typical Working Concentration:	27.22 g/L
Storage Temp:	2-6° C
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	Same formulation as B141 with 0.2 mg/L BA Unadjusted pH = 5.0 – 6.0 BM-2 differs from BM Medium in that it (BM-2) has no inorganic nitrogen and it contains hydrolyzed casein and L-glutamine, and BA.

Formula (mg/L)

Boric Acid	10
Cobalt Chloride•6H ₂ O	0.025
Cupric Sulfate•5H ₂ O	0.025
Na ₂ EDTA•2H ₂ O	37.25
Ferrous Sulfate•7H ₂ O	27.85
Magnesium Sulfate, Anhydrous	100
Manganese Sulfate•H ₂ O	25
Molybdc Acid (Sodium Salt) •2H ₂ O	0.25
Potassium Phosphate, Monobasic	300
Zinc Sulfate•7H ₂ O	10
6-Benzylaminopurine (BA)	0.2

D-Biotin	0.05
Casein, Enzymatic Hydrolysate	500
Folic Acid	0.5
L-Glutamine	100
Glycine (Free Base)	2
myo-Inositol	100
Nicotinic Acid (Free Acid)	5
Pyridoxine•HCl	0.5
Sucrose	20,000
Thiamine•HCl	0.5

Application Notes

Plant Tissue Culture Tested

Well suited for the terrestrial orchids *Paphiopedilum* and *Phragmipedium*.

References

According to "The Orchid Seedbank Project" this was originally cited by Van Waes, 1984. In vitro studie van de kiemingsfysiologie van Westeuropese orchideeën. Thesis. Rijkuniversiteit Gent.

Revised 3/2012