

NEW



A great new tool for Gene Cloning

GelX: Rapid, precise and contamination free excision of spots & bands from gels

Newly available from Cleaver Scientific, GelX disposable gel excision tips. No more scalpels required, these tips offer a unique and convenient way to excise bands from gels. GelX tips allow complete single handed and contamination free removal of the gel band of interest. The gel band is neatly extracted and contained in the tip. It can then be simply dispensed into a tube and the used tip ejected and disposed of.

- Eliminates scalpel damage to transilluminator or gel tray
- Designed for removing bands from agarose and protein gels
- One handed operation
- Push button gel release
- Push button tip release
- Avoid cross contamination
- Easy gel slice transfer GelX tips fit standard 1000 μ l pipettors
- Autoclavable
- DNase and RNase free
- Sterilized versions are also available



Ordering information:-

CSL-gelX4
4mm x 1mm
250/ bag

CSL- gelX 4RACK
4mm x 1mm
5 racks of 48

CSL- gelX 6.5
6.5mm x 1mm
250/ bag

CSL- gelX6.5RACK
6.5mm x 1mm
5 racks of 48