



**ELISA Coating Buffer  
E107  
Carbonate-Bicarbonate Buffer  
Capsules**

**Catalog No.**      **Size**  
**E107**                **2 capsules**

**Description:**

Carbonate-bicarbonate buffer capsules

Sigma Chemical #C3041

Package of 2

**Application:**

Capsules are suitable for use in protein coating procedures on multiwell plates or plastic tubes.

**Preparation:**

Dissolve contents of one capsule in 100 ml deionized water to yield a 0.05 M carbonate-bicarbonate buffer, pH 9.6, at 25 °C

Do not add the gelatin capsule to your solution.

**Storage:**

Store prepared coating buffer at room temperature for up to one year.

For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals.

Not for animal or human consumption.

**Related Products**

TMB Substrate	E102
Postcoat Solution	E104
Microtiter plates	E105
Wash Solution	E106
ELISA Accessory Starter Kit I	E101
ELISA Accessory Starter Kit II	E103