

Normal Porcine Blocking Buffer
(10% Normal Porcine Serum)

Storage: 2-8°C

Catalog No: E32-18 18 ml
E32-100 100ml

Intended Use:

Normal Porcine Blocking buffer is 10% normal porcine serum solution. Normal porcine serum is used to block non-specific staining when a porcine second antibody is used during immunohistochemical staining procedure.

Kit Components:

Catalog No.	Size
E32-18	18 ml (ready to use)
E32-100	100ml (ready to use)

Solution contains 0.05% NaN₃.

Recommended Protocol:

1. Prior to application of primary antibody apply Normal Porcine Blocking buffer on tissue slides and incubate the tissue section for 10 minutes at room temperature
2. Blot off Normal Porcine Blocking buffer (do not rinse) and add primary antibody

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

Storage:

Store at 2-8°C.

Remarks:

For research use only.