

Data sheet



Cat. No.	Amount
PR-BA132-S	250 μ g
PR-BA132-L	1 mg

For *in vitro* use only.

Quality guaranteed for 12 months.

Store at -80°C.

Avoid repeated freeze-thaw cycles.

Culture system:

Liquid culture

Protein concentration:

Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.

Preservation:

None

Analysis:

ELISA

Preparation:

Disrupted somatic material is clarified by centrifugation in phosphate buffered saline.

Inactivation:

Viable organisms are removed by centrifugation. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.