

Sterilising autoclaves

CertoClav



Aluminium alloy, not acid resistant

The compact CertoClav can be used in even the smallest laboratory. Everything that is constantly needed can therefore be sterilised at the workplace; cultures, glass containers, metal, rubber, textile and steam resistant synthetic materials.

The natural flow of steam generated by the round chamber, and the design of the steam release cock, guarantee effective venting, the prerequisite of perfect sterilisation of hollow and cylindrical shaped items.

CV-EL models: With four operating temperatures/pressures: 115 °C/0,7 bar; 121 °C/1,3 bar; 125 °C/1,4 bar; 140 °C/2,7 bar

Multicontrol models: Microprocessor operated and freely programmable; two fixed and seven free programmes. Fractionated venting. Separate control unit with digital display and RS232 serial interface.

Model	CV-EL 12 L GS	CV-EL 18 L GS	Multicontrol 12 L	Multicontrol 18 L
Capacity (l)	12	18	12	18
Temperature (°C)	140			
Heating power (W)	1,9			
Max. pressure (bar)	2,7			
Internal Ø×H (mm)	240×260	240×400	240×260	240×400
W×D×H of control unit (mm)	-		240×160×90	
Weight (kg)	6,1	7,7		9,3
Weight of control unit (kg)	-		1,6	

Supplied with instrument tray and tripod in aluminium.

Cena netto: 12922 zł

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium